



Body Fluid Identification Test No. 11-579 Summary Report

This proficiency test was sent to 96 participants. Each participant received a sample pack consisting of two known bloodstains and two questioned stains which they were requested to analyze using their existing protocols. Data were returned from 75 participants (78% response rate) and are compiled into the following tables:

	<u>Page</u>
<u>Manufacturer's Information</u>	<u>2</u>
<u>Summary Comments</u>	<u>3</u>
<u>Key to Screening Results</u>	<u>4</u>
<u>Table 1: Serology Screening Results</u>	<u>5</u>
<u>Table 2: Screening Conclusions</u>	<u>93</u>
<u>Table 3: Additional Comments</u>	<u>117</u>
<u>Appendix: Data Sheet</u>	<u>121</u>

This report contains the data received from the participants in this test. Since these participants are located in many countries around the world, and it is their option how the samples are to be used (e.g., training exercise, known or blind proficiency testing, research and development of new techniques, etc.), the results compiled in the Summary Report are not intended to be an overview of the quality of work performed in the profession and cannot be interpreted as such. The Summary Comments are included for the benefit of participants to assist with maintaining or enhancing the quality of their results. These comments are not intended to reflect the general state of the art within the profession.

Participant results are reported using a randomly assigned "WebCode". This code maintains participant's anonymity, provides linking of the various report sections, and will change with every report.

Manufacturer's Information

Each sample pack consisted of six items for testing, four of which were swabs and two of which were clothing items. Participants were asked to examine Items 1 and 2 for the body fluids that they typically encounter in casework. Items 3-6 were to be examined for scenario specific body fluids.

Each sample pack contained a package of unused swabs, questioned swabs (Items 1, 2, 3 and 4) and two questioned items of clothing, a glove (Item 5) and a white tank top (Item 6). The stain on Item 1 was prepared using pooled human saliva. Items 2, 4, 5 and 6 were prepared from human whole blood which was drawn into citric acid preservative bloodbank bags. The stain on Item 3 contained semen. The semen, procured from a cryobank, was from a single donor, but stored frozen in multiple vials at -196° C in liquid nitrogen tanks. The semen was thawed and mixed 1:1 with TAE buffer. This mixture was used to spot the Item 3 swabs. Stains from different sources were prepared at separate times and were packaged once they were thoroughly dried.

Sample Preparation Schedule

<u>Item</u>	<u>Substrate</u>	<u>Fluid</u>	<u>Volume</u>
1	4 Swabs	Saliva	100 µl/swab
2	4 Swabs	Blood	100 µl/swab
3	4 Swabs	Semen	100 µl/swab
4	4 Swabs	Blood	100 µl/swab
5	Glove	Blood	100 µl
6	White Tank Top	Blood	100 µl

Completed sample sets were stored at -20°C until shipment on September 7, 2011, following receipt of predistribution results.

Summary Comments

The 11-579 Body Fluid Identification test was designed to allow participants to assess their proficiency in the identification of stains that may contain dried body fluids. Each participant received six items for testing, four of which were packages of swabs and two were items of clothing. No scenario was provided for Items 1 and 2; participants were requested to test for body fluids that they typically encounter in casework. Scenarios were provided for the remaining items; participants were asked to treat these as they would in casework.

Item 1 consisted of four swabs containing saliva. Of the participants that reported a conclusion for saliva, the majority of these were positive, exclusively using presumptive testing. One participant reported negative results using a presumptive starch iodine test and one participant was inconclusive using a presumptive amylase test. Most participants reported negative results for both blood and semen testing, however eight reported "positive" or "indicated" conclusions for blood based upon positive results with the Hematrace test.

Item 2 consisted of four swabs containing blood. Of the participants that reported a conclusion for blood all were positive. The majority of these participants used a combination of presumptive and confirmatory tests. All participants reported negative conclusions for semen and most also for saliva, however one concluded positive for saliva based on a positive RSID test and five concluded saliva was possible based on starch/iodine test results.

Item 3 consisted of four swabs containing semen. Of the participants that reported a conclusion for semen all were positive. The majority of these participants used a combination of presumptive and confirmatory tests. Most participants reported negative results for both blood and saliva testing, however one concluded positive for blood based on a positive Hematrace test and four were positive for saliva based on positive RSID or "salivary amylase" test results.

Item 4 consisted of four swabs containing blood. Of the participants reported a conclusion for blood all were positive. The majority of these participants used a combination of presumptive and confirmatory tests. Most participants reported negative results for both semen and saliva testing, however three concluded either "positive" or "indicated" based on positive PSA test results and one concluded positive for saliva based on a positive RSID test.

Item 5 consisted of one glove containing blood. Of the participants that reported a conclusion for blood all were positive. The majority of these participants used a combination of presumptive and confirmatory tests. Of the participants that reported conclusions for semen all were negative and all but one were also negative for saliva. The one who was positive for saliva based it on a positive RSID test.

Item 6 consisted of one white tank top containing blood. Of the participants that reported a conclusion for blood all were positive. The majority of these participants used a combination of presumptive and confirmatory tests. One participant performed three presumptive tests and two of these were negative for blood but their overall conclusion was positive. Most of the participants that reported conclusions for semen and saliva were negative, however one reported positive for semen based on a positive ALS result and another reported positive for saliva based on a positive RSID test.

Key for Screening Tests Used

Participants were asked to use, where possible, the following chart of abbreviated screening test names. This was not an all inclusive list and was not designed to determine what tests should be performed. Participants were advised that tests not on this list may be used for screening.

Test	Abbreviation
Acid Phosphatase	AP
Kastle Meyer	KM
Leucomalachite Green	LMG
Microscopic	Micro
Ortho-tolidine	O-tol
Phenolphthalein-Tetramethyl benzidin	PTMB
Prostate Specific Antigen	PSA
Tetramethyl benzidine	TMB

Serology Screening Results

Indicate the results of any screening tests performed on the each Item.

TABLE 1 - Item 1

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
2XTJD9	Hemastix	Negative	Presumptive
	O-tol	Negative	Presumptive
4F8HWC	TMB	Negative	Presumptive
6DCQBC	KM	Negative	Presumptive
6RBUDH	KM	Negative	Presumptive
7A8K9C	KM	Negative	Presumptive
7N9F6H	LMG	Negative	Presumptive
	Hematrace	Positive	Confirmatory
7T4D8A	KM	Negative	Presumptive
7Y6RN4	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
8D22XB	KM	Negative	Presumptive
8WXUW9	Visual	Negative	Presumptive
927828	KM	Negative	Presumptive
AWM2QV	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
BLJ79Q	O-tol	Negative	Presumptive
BUXV7Y	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
C74XV9	KM	Negative	Presumptive
CT4Y9C	KM	Negative	Presumptive
CUWH34	KM	Negative	Presumptive
CWLNMZ	O-tol	Negative	Presumptive
	RSID (Rapid Stain Identification Kit)	Negative	Confirmatory
CY79E6	KM	Negative	Presumptive

TABLE 1 - Item 1

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
D8Y9K2	KM	Negative	Presumptive
DBGZH6	KM	Negative	Presumptive
DBYGLQ	KM	Negative	Presumptive
EBXDL8	KM	Negative	Presumptive
ELA9PQ	Visual	Negative	Presumptive
EXNXE3	KM	Negative	Presumptive
F92HC4	kastle mayer[sic]	Negative	Presumptive
FB83Y4	o-tol	Negative	Presumptive
FFJD6X	KM	Negative	Presumptive
FGBE6U	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
FPTPHV	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
GE84NK	Ortho-tolidine (O-tol)	Negative	Presumptive
GUWFC9	LMG	Negative	Presumptive
	Hematrace	Positive	Confirmatory
GZ76ZN	Visual	Negative	Presumptive
H4ZKQM	Visual	Negative	Presumptive
H7PRBJ	o-Tol	Negative	Presumptive
HAPXUW	KM	Negative	Presumptive
HJCUHX	KM	Negative	Presumptive
J8DD62	KM	Negative	Presumptive
JVLQY7	KM	Negative	Presumptive
K28BTQ	KM	Negative	Presumptive
	ABACard Hematrace	Negative	Confirmatory
K3WX4T	KM	Negative	Presumptive
KEA3BP	KM	Negative	Presumptive
	ABACard Hematrace	Negative	Confirmatory

TABLE 1 - Item 1

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
KETDZ6	LMG	Negative	Presumptive
	Hematrace	Positive	Confirmatory
KGHKK2	KM	Negative	Presumptive
KJN692	KM	Negative	Presumptive
	Characteristic Appearance	Negative	Presumptive
L9N8AL	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
LQLPYG	O-tol	Negative	Presumptive
LRV4H3	LMG	Negative	Presumptive
	Hematrace	Positive	Confirmatory
LUKB3Z	KM	Negative	Presumptive
M3B3FN	TMB (hemastix)	Positive	Presumptive
	Hematrace	Positive	Confirmatory
M48Q7U	KM	Negative	Presumptive
MGAGMT	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
NALRPG	KM	Negative	Presumptive
PF8YRL	Tetramethylbenzidine (TMB)	Negative	Presumptive
R7EA6H	KM	Negative	Presumptive
	ABAcad Hematrace	Negative	Confirmatory
R7YFFN	KM	Negative	Presumptive
RJZCCV	LMG	Negative	Presumptive
	HemaTrace	Positive	Confirmatory
T9XLJE	Visual	Negative	
TGBTYH	O-tol	Negative	Presumptive
	Hematrace	Negative	Confirmatory
TGTPMB	KM	Negative	Presumptive
ULKDJ9	o-Tol	Negative	Presumptive
VK47JR	KM	Negative	Presumptive

TABLE 1 - Item 1

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
WAZNX9	Tetrabase test	Negative	Confirmatory
	RSID blood	Negative	Presumptive
WQ4H49	Visual	Negative	Presumptive
WQNYJB	KM	Negative	Presumptive
XT8BX7	KM	Negative	Presumptive
XXLHBR	LMG	Negative	Presumptive
	Hematrace	Positive	Confirmatory
Y7WPEE	KM	Negative	Presumptive
ZG8GHE	KM	Negative	Presumptive
ZHYHGB	Hematrace	Negative	Confirmatory
	KM	Negative	Presumptive
ZVKJNN	LMG	Negative	Presumptive
	Hematrace	Positive	Confirmatory
ZW33KZ	O-tol	Negative	Presumptive

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
KM	46	Negative
Hematrace	21	Negative

TABLE 1 - Item 1

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
2XTJD9	ALS - Polilight	Negative	Presumptive
	ALS - Crime-lite	Negative	Presumptive
4F8HWC	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
6DCQBC	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
6RBUDH	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
7A8K9C	AP	Inconclusive	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
7N9F6H	AP	Negative	Presumptive
	micro	Negative	Confirmatory
	P30 (for PSA)	Negative	Confirmatory
7T4D8A	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
7Y6RN4	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
8D22XB	AP	Inconclusive	Presumptive
	P30 PSA	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
8WXUW9	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
927828	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory

TABLE 1 - Item 1

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
AWM2QV	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
BLJ79Q	AP	Positive	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory
BUXV7Y	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
C74XV9	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
CT4Y9C	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro (KPIC)	Negative	Confirmatory
CUWH34	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
CWLNMZ	AP	Negative	Presumptive
	RSID	Positive	Confirmatory
CY79E6	AP	Negative	Presumptive
	PSA (ABAcad, P30)	Negative	Presumptive
	Micro (KPIC)	Negative	Confirmatory
D8Y9K2	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
DBGZH6	AP	Negative	Presumptive
	PSA (p30)	Negative	Presumptive
	Microscopic	Negative	Confirmatory
DBYGLQ	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
EBXDL8	Acid Phosphatase (AP)	Negative	Presumptive
ELA9PQ	AP	Negative	Presumptive

TABLE 1 - Item 1

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
EXNXE3	AP	Inconclusive	Presumptive	
	PSA (ABACard)	Negative	Presumptive	
	Micro	Negative	Confirmatory	
F92HC4	Acid Phosphatase	Negative	Presumptive	
	Choline/Florence Iodine	Negative	Presumptive	
FB83Y4	Acid phosphatase (AP)	Negative	Presumptive	
	microscopy (sperm)	Negative	Confirmatory	
FFJD6X	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
FGBE6U	AP	Negative	Presumptive	
	ABACard P30	Negative	Presumptive	
	micro	Negative	Confirmatory	
FPTPHV	AP	Negative	Presumptive	
	ABACard P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
GE84NK	Alternate Light Source (ALS)	Positive	Presumptive	
	Acid Phosphatase (AP)	Negative	Presumptive	
GUWFC9	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
GZ76ZN	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	Rsid Semen	Negative	Presumptive	
H4ZKQM	AP	Negative	Presumptive	
H7PRBJ	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
HAPXUW	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	

TABLE 1 - Item 1

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
HJCUHX	AP	Negative	Presumptive	
	P30/PSA	Negative	Presumptive	
	KPIC/Micro	Negative	Confirmatory	
J8DD62	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
JVLQY7	AP	Negative	Presumptive	
K28BTQ	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
K3WX4T	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
KEA3BP	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
KETDZ6	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	P30	Negative	Confirmatory	
KGHKK2	AP	Negative	Presumptive	
	PSA	Negative	Confirmatory	
KJN692	AP	Negative	Presumptive	
	ALS	Positive	Presumptive	
L9N8AL	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
LQLPYG	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
LRV4H3	AP	Negative	Presumptive	
	Microscopic	Negative	Confirmatory	
	PSA	Negative	Confirmatory	

TABLE 1 - Item 1

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
LUKB3Z	AP	Negative	Presumptive
	Microscopic examination	Negative	Confirmatory
M3B3FN	Micro	Negative	Confirmatory
	AP	Negative	Presumptive
M48Q7U	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro - KPIC	Negative	Confirmatory
MGAGMT	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
	Micro (CT)	Negative	Confirmatory
N4Y7VW	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
NALRPG	AP	Negative	Presumptive
PF8YRL	Microscopy	Negative	Confirmatory
	Acid Phosphatase (AP)	Negative	Presumptive
R7EA6H	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
R7YFFN	AP	Negative	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
RJZCCV	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory
T9XLJE	AP	Negative	Presumptive
TGBTYH	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
TGTPMB	AP	Negative	Presumptive
	Micro	Negative	Confirmatory

TABLE 1 - Item 1

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
ULKDJ9	ALS	Positive	Presumptive	
	AP	Negative	Presumptive	
	ABAcad P30	Negative	Confirmatory	
	Sperm Cell Extraction/Micro	Negative	Confirmatory	
VK47JR	AP	Inconclusive	Presumptive	
	PSA	Negative	Presumptive	
	Micro	Negative	Confirmatory	
WAZNX9	Microscopic analysis	Negative	Confirmatory	
	PSA test	Negative	Presumptive	
	RSID semen	Negative	Presumptive	
WLAQ8R	ALS (Crime lite)	Negative	Presumptive	
	Phosphatesmo KM	Negative	Confirmatory	
WQ4H49	Acid Phosphatase	Negative	Presumptive	
WQNYJB	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
XT8BX7	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
XXLHBR	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
Y7WPEE	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
ZG8GHE	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
ZHYHGB	ABAcad P30	Negative	Presumptive	
	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
ZVKJNN	AP	Negative	Presumptive	
	Microslide	Negative	Confirmatory	
	PSA	Negative	Confirmatory	

TABLE 1 - Item 1

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
ZW33KZ	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
AP	71	Negative
Microscopy	58	Negative

TABLE 1 - Item 1

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	Phadebas	Positive	Presumptive
4F8HWC	Phadebas amylase test	Positive	Presumptive
	Micro	Positive	Presumptive
6DCQBC	Phadebas	Positive	Presumptive
6RBUDH	Phadebas	Positive	Presumptive
7A8K9C	Phadebas	Positive	Presumptive
7N9F6H	Starch/Iodine	Positive	Presumptive
7T4D8A	Phadebas Press Test	Positive	Presumptive
7Y6RN4	Phadebas	Positive	Presumptive
8D22XB	Phadebas	Positive	Presumptive
8WXUW9	RSID	Positive	Confirmatory
927828	Phadebas	Positive	Presumptive
AWM2QV	Phadebas	Positive	Presumptive
BLJ79Q	RSID	Positive	Confirmatory
BUXV7Y	Phadebas	Positive	Presumptive
C74XV9	Phadebas	Positive	Presumptive
CT4Y9C	Phadebas	Positive	Presumptive
CUWH34	Phadebas press test	Positive	Presumptive
CWLNMZ	RSID	Positive	Confirmatory
CY79E6	Phadebas	Positive	Presumptive
D8Y9K2	Phadebas	Positive	Presumptive
DBGZH6	Phadebas	Positive	Presumptive
DBYGLQ	Amylase Radial Diffusion	Positive	Presumptive
EBXDL8	Phadebas Tube test	Positive	Presumptive
EXNXE3	Phadebas	Positive	Presumptive
F92HC4	Phadebas	Positive	Presumptive
FB83Y4	amylase diffusion	Positive	Presumptive
	microscopy (NEC)	Positive	Presumptive

TABLE 1 - Item 1

Saliva				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
FFJD6X	Phadebas	Positive	Presumptive	
FGBE6U	Phadebas	Positive	Presumptive	
FPTPHV	Phadebas	Positive	Presumptive	
GE84NK	RSID - Saliva	Positive	Confirmatory	
GUWFC9	Starch-Iodine	Positive	Presumptive	
GZ76ZN	Rsid Saliva	Positive	Presumptive	
H7PRBJ	RSID	Positive	Confirmatory	
HAPXUW	Phadebas	Positive	Presumptive	
HJCUHX	Phadebas	Positive	Presumptive	
J8DD62	Phadebas	Positive	Presumptive	
JVLQY7	Phadebas	Positive	Presumptive	
K28BTQ	Phadebas	Positive	Presumptive	
K3WX4T	Phadebas	Positive	Presumptive	
KEA3BP	Phadebas	Positive	Presumptive	
KETDZ6	Starch Iodine	Positive	Presumptive	
KGHKK2	a Amylase	Positive	Confirmatory	
KJN692	ALS	Positive	Presumptive	
	RSID - Saliva	Positive	Confirmatory	
L9N8AL	Phadebas	Positive	Presumptive	
LQLPYG	RSID-Saliva	Positive	Confirmatory	
LRV4H3	Starch-Iodine	Positive	Presumptive	
LUKB3Z	Phadebas	Positive	Presumptive	
M3B3FN	RSID - saliva	Positive	Confirmatory	
M48Q7U	Phadebas	Positive	Presumptive	
MGAGMT	Phadebas	Positive	Presumptive	
	Amilase[sic] - RSID	Positive	Confirmatory	
N4Y7WW	Amylase Gel	Positive	Presumptive	
	Micro	Positive	Presumptive	

TABLE 1 - Item 1

Saliva				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
NALRPG	Phadebas Amylase	Positive	Presumptive	
PF8YRL	Phadebas	Positive	Presumptive	
R7EA6H	Phadebas	Positive	Presumptive	
R7YFFN	Phadebas	Positive	Presumptive	
RJZCCV	Starch-Iodine	Positive	Presumptive	
TGBTYH	Salivary amylase	Positive	Presumptive	
TGTPMB	Radial Diffusion	Positive	Presumptive	
ULKDJ9	ALS	Positive	Presumptive	
	RSID-Saliva	Positive	Confirmatory	
VK47JR	Phadebas	Positive	Presumptive	
WAZNX9	RSID saliva	Positive	Presumptive	
WLAQ8R	RSID Saliva	Positive	Confirmatory	
WQNYJB	Phadebas	Positive	Presumptive	
XT8BX7	Amylase	Inconclusive	Presumptive	
XXLHBR	Starch Iodine Test	Negative	Presumptive	
Y7WPEE	Phadebas	Positive	Presumptive	
ZG8GHE	Phadebas	Positive	Presumptive	
ZHYHGB	Phadebas	Positive	Presumptive	
ZVKJNN	Starch Iodine	Positive	Presumptive	
ZW33KZ	RSID-saliva	Positive	Confirmatory	

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
Phadebas	43	Positive
RSID-Saliva	14	Positive

TABLE 1 - Item 1

Other Body Fluid Types			
Webcode	Body Fluid Type - Test Performed	Results	Presumptive/Confirmatory
7N9F6H	epithelial cells - possible; swabbed for poss. e.c.'s		Presumptive
GUWFC9	Epithelial cells - Micro	Positive	Confirmatory
GZ76ZN	Urine - Rsid Urine	Negative	
KGHKK2	Semen - ALS	Negative	Presumptive
LRV4H3	Epithelial cells - Microscopic	Positive	Confirmatory
MGAGMT	ALS	Positive	Presumptive
	Cells - Micro	Positive	Confirmatory
NALRPG	Urine - Creatinine	Negative	Presumptive
	Urine - DMAC	Negative	Presumptive
RJZCCV	Epithelial cells - Micro	Positive	Confirmatory
WLAQ8R	Urine - RSID Urine	Negative	Confirmatory
XXLHBR	Epithelial cells - Micro.	Positive	Confirmatory
ZVKJNN	Epithelial Cell - Microslide	Positive	Confirmatory

TABLE 1 - Item 2

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
2XTJD9	Hemastix	Positive	Presumptive
	Otol	Positive	Presumptive
4F8HWC	TMB	Positive	Presumptive
	Haematrace[sic]	Positive	Confirmatory
6DCQBC	KM	Positive	Presumptive
6RBUDH	KM	Positive	Presumptive
	Takayama	Positive	Confirmatory
	Ouchterlony-human	Positive	Confirmatory
7A8K9C	KM	Positive	Presumptive
7N9F6H	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7T4D8A	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7Y6RN4	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
8D22XB	KM	Positive	Presumptive
8WXUW9	visual	Positive	Presumptive
	Otol	Positive	Presumptive
927828	KM	Positive	Presumptive
AWM2QV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BLJ79Q	O-tol	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BUXV7Y	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
C74XV9	KM	Positive	Presumptive
CT4Y9C	KM	Positive	Presumptive

TABLE 1 - Item 2

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
CUWH34	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
CWLNMZ	O-tol	Positive	Presumptive
	RSID	Positive	Confirmatory
CY79E6	KM	Positive	Presumptive
D8Y9K2	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
DBGZH6	KM	Positive	Presumptive
DBYGLQ	KM	Positive	Presumptive
	Hematrace	Positive	Presumptive
EBXDL8	KM	Positive	Presumptive
ELA9PQ	Visual	Positive	Presumptive
	LMG	Positive	Presumptive
	ABA Hema Trace	Positive	Confirmatory
EXNXE3	KM	Positive	Presumptive
F92HC4	kastle mayer[sic]	Positive	Presumptive
FB83Y4	o-tol	Positive	Presumptive
	hematrace	Positive	Presumptive
FFJD6X	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
FGBE6U	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
FPTPHV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
GE84NK	O-tol	Positive	Presumptive
	HemaTrace (HT)	Positive	Confirmatory
GUWFC9	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory

TABLE 1 - Item 2

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
GZ76ZN	Visual	Positive	Presumptive
	Km	Positive	Presumptive
	Haematrace[sic]	Positive	Presumptive
H4ZKQM	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
	Visual	Positive	Presumptive
H7PRBJ	o-Tol	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
HAPXUW	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
HJCUHX	KM	Positive	Presumptive
J8DD62	KM	Positive	Presumptive
JVLQY7	KM	Positive	Presumptive
K28BTQ	KM	Positive	Presumptive
	ABAcad Hematrace	Positive	Confirmatory
K3WX4T	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
KEA3BP	KM	Positive	Presumptive
	ABAcad Hematrace	Positive	Confirmatory
KETDZ6	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
KGHKK2	KM	Positive	Presumptive
KJN692	KM	Positive	Presumptive
	Characteristic Appearance	Positive	Presumptive
L9N8AL	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
LQLPYG	O-tol	Positive	Presumptive
	HemaTrace	Positive	Confirmatory

TABLE 1 - Item 2

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
LRV4H3	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
LUKB3Z	KM	Positive	Presumptive	
M3B3FN	TMB (hemastix)	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
M48Q7U	KM	Positive	Presumptive	
MGAGMT	KM	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
N4Y7VW	O-Tol	Positive	Presumptive	
	Hematrace	Positive	Presumptive	
NALRPG	KM	Positive	Presumptive	
	Micro (Mod Tak)	Positive	Confirmatory	
PF8YRL	TMB	Positive	Presumptive	
R7EA6H	KM	Positive	Presumptive	
	ABAcad Hematrace	Positive	Confirmatory	
R7YFFN	KM	Positive	Presumptive	
RJZCCV	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
T9XLJE	Hemochromogen	Positive	Confirmatory	
TGBTYH	O-tol	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
TGTPMB	KM	Positive	Presumptive	
	Hematrace	Positive	Presumptive	
ULKDJ9	o-Tol	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
VK47JR	KM	Positive	Presumptive	
WAZNX9	Tetrabase test	Positive	Confirmatory	
	RSID blood	Positive	Presumptive	

TABLE 1 - Item 2

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
WLAQ8R	Hemastix	Positive	Presumptive	
	Hexagon OBTI	Positive	Confirmatory	
WQ4H49	Visual	Positive	Presumptive	
	LMG	Positive	Presumptive	
	ABAcad Hematrace	Positive	Confirmatory	
WQNYJB	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
XT8BX7	KM	Positive	Presumptive	
	Hematrace	Positive	Presumptive	
XXLHBR	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
Y7WPEE	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
ZG8GHE	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
ZHYHGB	Hematrace	Positive	Confirmatory	
	KM	Positive	Presumptive	
ZVKJNN	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
ZW33KZ	O-tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
KM	47	Positive
Hematrace	46	Positive

TABLE 1 - Item 2

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
2XTJD9	ALS - Polilight	Negative	Presumptive
	ALS - Crime lite	Negative	Presumptive
4F8HWC	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
6DCQBC	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
6RBUDH	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
7A8K9C	AP	Negative	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
7N9F6H	AP	Negative	Presumptive
	micro	Negative	Confirmatory
	P30 (for PSA)	Negative	Confirmatory
7T4D8A	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
7Y6RN4	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
8D22XB	AP	Negative	Presumptive
	P30 PSA	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
8WXUW9	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
927828	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory

TABLE 1 - Item 2

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
AWM2QV	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
BLJ79Q	AP	Inconclusive	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory
BUXV7Y	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
C74XV9	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
CT4Y9C	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro (KPIC)	Negative	Confirmatory
CUWH34	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
CWLNMZ	AP	Negative	Presumptive
	RSID	Negative	Confirmatory
CY79E6	AP	Negative	Presumptive
	PSA (ABAcad, P30)	Negative	Presumptive
	Micro (KPIC)	Negative	Confirmatory
D8Y9K2	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
DBGZH6	AP	Inconclusive	Presumptive
	PSA (p30)	Negative	Presumptive
	Microscopic	Negative	Confirmatory
DBYGLQ	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
EBXDL8	Acid Phosphatase (AP)	Negative	Presumptive

TABLE 1 - Item 2

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
ELA9PQ	AP	Negative	Presumptive	
	ABA P30	Negative	Confirmatory	
	microscope	Negative	Confirmatory	
EXNXE3	AP	Inconclusive	Presumptive	
	PSA	Negative	Presumptive	
	Micro	Negative	Confirmatory	
F92HC4	Acid Phosphatase	Negative	Presumptive	
	Choline/Florence Iodine	Negative	Presumptive	
FB83Y4	Acid phosphatase	Negative	Presumptive	
	micro (sperm)	Negative	Confirmatory	
FFJD6X	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
FGBE6U	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
FPTPHV	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Microscopic	Negative	Confirmatory	
GE84NK	ALS	Negative	Presumptive	
	AP	Negative	Presumptive	
GUWFC9	AP	Inconclusive	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
GZ76ZN	Ap	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	Rsid Semen	Negative	Presumptive	
H4ZKQM	AP	Negative	Presumptive	
	PSA	Negative	Confirmatory	
	Micro	Negative	Confirmatory	

TABLE 1 - Item 2

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
H7PRBJ	AP	Inconclusive	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
HAPXUW	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
HJCUHX	AP	Negative	Presumptive	
	P30/PSA	Negative	Presumptive	
	KPIC/Micro	Negative	Confirmatory	
J8DD62	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
JVLQY7	AP	Negative	Presumptive	
K28BTQ	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
K3WX4T	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
KEA3BP	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
KETDZ6	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	P30	Negative	Confirmatory	
KGHKK2	AP	Negative	Presumptive	
	PSA	Negative	Confirmatory	
KJN692	ALS	Negative	Presumptive	
	AP	Negative	Presumptive	
L9N8AL	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	

TABLE 1 - Item 2

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
LQLPYG	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
LRV4H3	AP	Negative	Presumptive	
	Microscopic	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
LUKB3Z	AP	Negative	Presumptive	
	Microscopic examination	Negative	Confirmatory	
M3B3FN	Micro	Negative	Confirmatory	
	AP	Negative	Presumptive	
M48Q7U	AP	Negative	Presumptive	
	PSA	Negative	Presumptive	
	Micro - KPIC	Negative	Confirmatory	
MGAGMT	AP	Negative	Presumptive	
	PSA	Negative	Confirmatory	
	Micro (CT)	Negative	Confirmatory	
N4Y7VW	PSA	Negative	Presumptive	
	Micro	Negative	Confirmatory	
NALRPG	AP	Negative	Presumptive	
PF8YRL	Microscopy	Negative	Confirmatory	
	AP	Negative	Presumptive	
R7EA6H	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
R7YFFN	AP	Negative	Presumptive	
	P30 (PSA)	Negative	Presumptive	
	KPIC (Micro)	Negative	Confirmatory	
RJZCCV	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	

TABLE 1 - Item 2

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
T9XLJE	AP	Negative	Presumptive	
TGBTYH	AP	Negative	Presumptive	
	PSA	Negative	Confirmatory	
TGTPMB	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
ULKDJ9	ALS	Negative	Presumptive	
	AP	Negative	Presumptive	
	ABAcad P30	Negative	Confirmatory	
	Sperm Cell Extraction/Micro	Negative	Confirmatory	
VK47JR	AP	Negative	Presumptive	
	PSA	Negative	Presumptive	
	Micro	Negative	Confirmatory	
WAZNX9	Microscopic analysis	Negative	Confirmatory	
	PSA test	Negative	Presumptive	
	RSID semen	Negative	Presumptive	
WLAQ8R	ALS (Crime lite)	Negative	Presumptive	
WQ4H49	Acid Phosphatase	Negative	Presumptive	
	ABAcad P30	Negative	Confirmatory	
	Microscopic (sperm)	Negative	Confirmatory	
WQNYJB	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
XT8BX7	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
XXLHBR	AP	Negative	Presumptive	
	Micro.	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
Y7WPEE	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
ZG8GHE	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	

TABLE 1 - Item 2

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
ZHYHGB	ABAcad P30	Negative	Presumptive
	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
ZVKJNN	AP	Negative	Presumptive
	Microslide	Negative	Confirmatory
	PSA	Negative	Confirmatory
ZW33KZ	AP	Inconclusive	Presumptive
	micro	Negative	Confirmatory
	PSA	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
AP	71	Negative
Microscopy	61	Negative

TABLE 1 - Item 2

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	Phadebas	Negative	Presumptive
4F8HWC	Phadebas amylase test	Negative	Presumptive
6DCQBC	Phadebas	Negative	Presumptive
6RBUDH	Phadebas	Negative	Presumptive
7A8K9C	Phadebas	Negative	Presumptive
7N9F6H	Starch/Iodine	Positive	Presumptive
7T4D8A	Phadebas Tube Test	Inconclusive	Presumptive
7Y6RN4	Phadebas	Negative	Presumptive
8D22XB	Phadebas	Negative	Presumptive
8WXUW9	RSID	Negative	Confirmatory
927828	Phadebas	Negative	Presumptive
AWM2QV	Phadebas	Negative	Presumptive
BLJ79Q	RSID	Negative	Confirmatory
BUXV7Y	Phadebas	Negative	Presumptive
C74XV9	Phadebas	Negative	Presumptive
CT4Y9C	Phadebas	Negative	Presumptive
CUWH34	Phadebas tube test	Inconclusive	Presumptive
CWLNMZ	RSID	Positive	Confirmatory
CY79E6	Phadebas	Negative	Presumptive
D8Y9K2	Phadebas	Negative	Presumptive
DBGZH6	Phadebas	Negative	Presumptive
DBYGLQ	Amylase Radial Diffusion	Negative	Presumptive
EBXDL8	Phadebas tube test	Negative	Presumptive
EXNXE3	Phadebas	Negative	Presumptive
F92HC4	Phadebas	Inconclusive	Presumptive
FB83Y4	amylase diffusion	Negative	Presumptive
	micro (NEC)	Negative	Presumptive
FFJD6X	Phadebas	Inconclusive	Presumptive

TABLE 1 - Item 2

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
FGBE6U	Phadebas	Negative	Presumptive
FPTPHV	Phadebas	Negative	Presumptive
GE84NK	RSID - Saliva	Negative	Confirmatory
GUWFC9	Starch-Iodine	Positive	Presumptive
GZ76ZN	Rsid Saliva	Negative	Presumptive
H7PRBJ	RSID	Negative	Confirmatory
HAPXUW	Phadebas	Negative	Presumptive
HJCUHX	Phadebas	Negative	Presumptive
J8DD62	Phadebas	Negative	Presumptive
JVLQY7	Phadebas	Negative	Presumptive
K28BTQ	Phadebas	Negative	Presumptive
K3WX4T	Phadebas	Negative	Presumptive
KEA3BP	Phadebas	Negative	Presumptive
KETDZ6	Starch Iodine	Positive	Presumptive
KGHKK2	α -Amylase	Negative	Confirmatory
KJN692	ALS	Negative	Presumptive
	RSID - Saliva	Negative	Confirmatory
L9N8AL	Phadebas	Negative	Presumptive
LQLPYG	RSID-Saliva	Negative	Confirmatory
LRV4H3	Starch-Iodine	Negative	Presumptive
LUKB3Z	Phadebas	Negative	Presumptive
M3B3FN	RSID - saliva	Negative	Confirmatory
M48Q7U	Phadebas	Negative	Presumptive
MGAGMT	Phadebas	Negative	Presumptive
	Amilase[sic] - RSID	Negative	Confirmatory
N4Y7WW	Amylase Gel	Negative	Presumptive
	Micros	Negative	Presumptive
NALRPG	Phadebas Amylase	Negative	Presumptive

TABLE 1 - Item 2

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
PF8YRL	Phadebas	Negative	Presumptive
R7EA6H	Phadebas	Negative	Presumptive
R7YFFN	Phadebas	Negative	Presumptive
RJZCCV	Starch-Iodine	Positive	Presumptive
TGBTYH	Salivary amylase	Negative	Presumptive
TGTPMB	Radial Diffusion	Negative	Presumptive
ULKDJ9	ALS	Negative	Presumptive
	RSID-Saliva	Negative	Confirmatory
VK47JR	Phadebas	Negative	Presumptive
WAZNX9	RSID saliva	Negative	Presumptive
WLAQ8R	RSID Saliva	Negative	Confirmatory
WQNYJB	Phadebas	Negative	Presumptive
XT8BX7	Amylase	Negative	Presumptive
XXLHBR	Starch Iodine	Positive	Presumptive
Y7WPEE	Phadebas	Inconclusive	Presumptive
ZG8GHE	Phadebas	Negative	Presumptive
ZHYHGB	Phadebas	Negative	Presumptive
ZVKJNN	Starch iodine	Negative	Presumptive
ZW33KZ	RSID-saliva	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
Phadebas	43	Negative
RSID-Saliva	14	Negative

TABLE 1 - Item 2

Other Body Fluid Types			
Webcode	Body Fluid Type - Test Performed	Results	Presumptive/Confirmatory
7N9F6H	epithelial cells - possible; swabbed for poss. e.c.'s		Presumptive
GUWFC9	Epithelial cells - Micro	Positive	Confirmatory
GZ76ZN	Urine - Rsid Urine	Negative	Presumptive
KGHKK2	Semen - ALS	Negative	Presumptive
LRV4H3	Epithelial cells - Microscopic	Positive	Confirmatory
MGAGMT	ALS	Negative	Presumptive
	Cells - Micro	Negative	Confirmatory
NALRPG	Urine - Creatinine	Negative	Presumptive
	Urine - DMAC	Negative	Presumptive
	Feces - Urobilinogen	Negative	Presumptive
	Species Origin (Animal)	Negative	Presumptive
RJZCCV	Epithelial Cells - Micro	Positive	Confirmatory
XXLHBR	Epithelial Cells - Micro.	Positive	Confirmatory
ZVKJNN	Epithelial Cell - Microslide	Positive	Confirmatory

TABLE 1 - Item 3

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
2XTJD9	Hemastix	Negative	Presumptive
	Otol	Negative	Presumptive
4F8HWC	TMB	Negative	Presumptive
6DCQBC	KM	Negative	Presumptive
6RBUDH	KM	Negative	Presumptive
7A8K9C	KM	Negative	Presumptive
7Y6RN4	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
8D22XB	KM	Negative	Presumptive
8WXUW9	visual	Negative	Presumptive
927828	KM	Negative	Presumptive
AWM2QV	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
BLJ79Q	O-tol	Negative	Presumptive
BUXV7Y	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
CT4Y9C	KM	Negative	Presumptive
CWLNMZ	O-tol	Negative	Presumptive
	RSID	Negative	Confirmatory
CY79E6	KM	Negative	Presumptive
D8Y9K2	KM	Negative	Presumptive
DBGZH6	KM	Negative	Presumptive
EBXDL8	KM	Negative	Presumptive
ELA9PQ	visual	Negative	Presumptive
EXNXE3	KM	Negative	Presumptive
F92HC4	kastle mayer[sic]	Negative	Presumptive

TABLE 1 - Item 3

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
FB83Y4	o-tol	Negative	Presumptive
FFJD6X	KM	Negative	Presumptive
FGBE6U	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
FPTPHV	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
GE84NK	O-tol	Negative	Presumptive
GZ76ZN	Visual	Negative	Presumptive
H4ZKQM	Visual	Negative	Presumptive
H7PRBJ	o-Tol	Negative	Presumptive
HAPXUW	KM	Negative	Presumptive
HJCUHX	KM	Negative	Presumptive
K28BTQ	KM	Negative	Presumptive
	ABAcad Hematrace	Negative	Confirmatory
K3WX4T	KM	Negative	Presumptive
KEA3BP	KM	Negative	Presumptive
	ABAcad Hematrace	Negative	Confirmatory
KGHKK2	KM	Negative	Presumptive
KJN692	Characteristic Appearance	Negative	Presumptive
	KM	Negative	Presumptive
L9N8AL	KM	Negative	Presumptive
	Hematrace	Negative	Confirmatory
LQLPYG	O-tol	Negative	Presumptive
M3B3FN	TMB (hemastix)	Negative	Presumptive
M48Q7U	KM	Negative	Presumptive
MGAGMT	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
NALRPG	KM	Negative	Presumptive

TABLE 1 - Item 3

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
PF8YRL	TMB	Negative	Presumptive
R7EA6H	KM	Negative	Presumptive
	ABACard Hematrace	Negative	Confirmatory
R7YFFN	KM	Negative	Presumptive
TGBTYH	O-tol	Negative	Presumptive
	Hematrace	Negative	Confirmatory
ULKDJ9	o-Tol	Negative	Presumptive
VK47JR	KM	Negative	Presumptive
WAZNX9	Tetrabase test	Negative	Confirmatory
	RSID blood	Negative	Presumptive
WQ4H49	Visual	Negative	Presumptive
WQNYJB	KM	Negative	Presumptive
Y7WPEE	KM	Negative	Presumptive
ZG8GHE	KM	Negative	Presumptive
ZHYHGB	Hematrace	Negative	Confirmatory
	KM	Negative	Presumptive
ZW33KZ	O-tol	Negative	Presumptive

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
KM	37	Negative
Hematrace	13	Negative

TABLE 1 - Item 3

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	AP	Positive	Presumptive
	ABAcad P30	Positive	Presumptive
	micro	Positive	Confirmatory
2XTJD9	ALS - Polilight	Positive	Presumptive
	ALS - Crimelite	Positive	Presumptive
4F8HWC	AP	Positive	Presumptive
	Micro	Positive	Confirmatory
6DCQBC	AP	Positive	Presumptive
	PSA	Positive	Presumptive
	Micro	Positive	Confirmatory
6RBUDH	AP	Positive	Presumptive
	Micro	Positive	Confirmatory
7A8K9C	AP	Positive	Presumptive
	P30 (PSA)	Positive	Presumptive
	KPIC (Micro)	Positive	Confirmatory
7N9F6H	AP	Positive	Presumptive
	micro	Positive	Confirmatory
	P30 (for PSA)	Positive	Confirmatory
7T4D8A	AP	Positive	Presumptive
	PSA	Positive	Confirmatory
7Y6RN4	AP	Positive	Presumptive
	ABAcad P30	Positive	Presumptive
	Micro	Positive	Confirmatory
8D22XB	AP	Positive	Presumptive
	P30 PSA	Positive	Presumptive
	KPIC (Micro)	Positive	Confirmatory
8WXUW9	AP	Positive	Presumptive
	Micro	Positive	Confirmatory

TABLE 1 - Item 3

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
927828	AP	Positive	Presumptive	
	PSA	Positive	Presumptive	
	Micro	Positive	Confirmatory	
AWM2QV	AP	Positive	Presumptive	
	ABAcad P30	Positive	Presumptive	
	Micro	Positive	Confirmatory	
BLJ79Q	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
	PSA	Positive	Confirmatory	
BUXV7Y	AP	Positive	Presumptive	
	ABAcad P30	Positive	Presumptive	
	Micro	Positive	Confirmatory	
C74XV9	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
CT4Y9C	AP	Positive	Presumptive	
	PSA	Positive	Presumptive	
	Micro (KPIC)	Positive	Confirmatory	
CUWH34	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
CWLNMZ	AP	Positive	Presumptive	
	Microscopy	Positive	Confirmatory	
	RSID	Positive	Confirmatory	
CY79E6	AP	Positive	Presumptive	
	PSA (ABAcad, P30)	Positive	Presumptive	
	Micro (KPIC)	Positive	Confirmatory	
D8Y9K2	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
DBGZH6	AP	Positive	Presumptive	
	PSA (p30)	Positive	Presumptive	
	Microscopic	Positive	Confirmatory	

TABLE 1 - Item 3

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
DBYGLQ	AP	Positive	Presumptive
	Micro	Positive	Confirmatory
EBXDL8	Acid Phosphatase (AP)	Positive	Presumptive
	Microscopic (Haematoxylin and Eosin stain)	Positive	Confirmatory
ELA9PQ	AP	Positive	Presumptive
	ABA P30	Positive	Confirmatory
	Microscope	Positive	Confirmatory
EXNXE3	AP	Positive	Presumptive
	PSA	Positive	Presumptive
	Micro	Positive	Confirmatory
F92HC4	Extraction of cells for spermatozoa	Positive	Confirmatory
	PSA	Positive	Presumptive
	Acid Phosphatase	Positive	Presumptive
FB83Y4	Acid phosphatase (AP)	Positive	Presumptive
	p30 (PSA)	Positive	Presumptive
	micro (sperm)	Positive	Confirmatory
FFJD6X	AP	Positive	Presumptive
	PSA	Positive	Confirmatory
FGBE6U	AP	Positive	Presumptive
	ABAcad P30	Positive	Presumptive
	Micro	Positive	Confirmatory
FPTPHV	AP	Positive	Presumptive
	ABAcad P30	Positive	Presumptive
	Micro	Positive	Confirmatory
GE84NK	ALS	Positive	Presumptive
	AP	Positive	Presumptive
	Sperm search (microscopic)	Positive	Confirmatory
GUWFC9	AP	Positive	Presumptive
	Micro	Positive	Confirmatory

TABLE 1 - Item 3

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
GZ76ZN	Ap	Positive	Presumptive	
	Micro	Positive	Confirmatory	
H4ZKQM	AP	Positive	Presumptive	
	PSA	Positive		
	Micro	Positive		
H7PRBJ	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
	PSA	Positive	Confirmatory	
HAPXUW	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
HJCUHX	AP	Positive	Presumptive	
	P30/PSA	Positive	Presumptive	
	KPIC/Micro	Positive	Confirmatory	
J8DD62	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
JVLQY7	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
K28BTQ	AP	Positive	Presumptive	
	ABAcad P30	Positive	Presumptive	
	Micro	Positive	Confirmatory	
K3WX4T	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
KEA3BP	AP	Positive	Presumptive	
	ABAcad P30	Positive	Presumptive	
	Micro	Positive	Confirmatory	
KETDZ6	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
KGHKK2	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
	Micro	Positive	Confirmatory	

TABLE 1 - Item 3

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
KJN692	ALS	Positive	Presumptive	
	AP	Positive	Presumptive	
	Microscopic (sperm)	Positive	Confirmatory	
L9N8AL	AP	Positive	Presumptive	
	ABAcad P30	Positive	Presumptive	
	micro	Positive	Confirmatory	
LQLPYG	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
	PSA	Positive	Confirmatory	
LRV4H3	AP	Positive	Presumptive	
	Microscopic	Positive	Confirmatory	
LUKB3Z	AP	Positive	Presumptive	
	Microscopic examination	Positive	Confirmatory	
M3B3FN	Micro	Positive	Confirmatory	
M48Q7U	AP	Positive	Presumptive	
	PSA	Positive	Presumptive	
	Micro - KPIC	Positive	Confirmatory	
MGAGMT	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
	Micro (CT)	Positive	Confirmatory	
N4Y7VW	PSA	Positive	Presumptive	
	Micro	Positive	Confirmatory	
NALRPG	AP	Positive	Presumptive	
	Micro (Xmas Tree Stain)	Positive	Confirmatory	
PF8YRL	Microscopy	Positive	Confirmatory	
R7EA6H	AP	Positive	Presumptive	
	ABAcad P30	Positive	Presumptive	
	Micro	Positive	Confirmatory	

TABLE 1 - Item 3

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
R7YFFN	AP	Positive	Presumptive	
	P30 (PSA)	Positive	Presumptive	
	KPIC (Micro)	Positive	Confirmatory	
RZCCV	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
T9XLJE	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
TGBTYH	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
TGTPMB	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
ULKDJ9	ALS	Positive	Presumptive	
	AP	Positive	Presumptive	
	ABAcad P30	Positive	Confirmatory	
	Sperm Cell Extraction/Micro	Positive	Confirmatory	
VK47JR	AP	Positive	Presumptive	
	PSA	Positive	Presumptive	
	Micro	Positive	Confirmatory	
WAZNX9	Microscopic analysis	Positive	Confirmatory	
WLAQ8R	ALS (Crime lite)	Positive	Presumptive	
	Phosphatesmo KM	Positive	Confirmatory	
WQ4H49	Acid Phosphatase	Positive	Presumptive	
	ABAcad P30	Positive	Confirmatory	
	Microscopic (sperm)	Positive	Confirmatory	
WQNYJB	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
XT8BX7	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	

TABLE 1 - Item 3

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
XXLHBR	AP	Positive	Presumptive	
	Micro.	Positive	Confirmatory	
	PSA	Positive	Confirmatory	
Y7WPEE	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
ZG8GHE	AP	Positive	Presumptive	
	PSA	Positive	Confirmatory	
ZHYHGB	ABACard P30	Positive	Presumptive	
	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	
ZVKJNN	AP	Positive	Presumptive	
	Microslide	Positive	Confirmatory	
ZW33KZ	AP	Positive	Presumptive	
	Micro	Positive	Confirmatory	

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
AP	69	Positive
Microscopy	62	Positive

TABLE 1 - Item 3

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	Phadebas	Negative	Presumptive
4F8HWC	Phadebas amylase test	Negative	Presumptive
6DCQBC	Phadebas	Negative	Presumptive
6RBUDH	Phadebas	Negative	Presumptive
7A8K9C	Phadebas	Negative	Presumptive
7T4D8A	Phadebas Press Test	Negative	Presumptive
7Y6RN4	Phadebas	Negative	Presumptive
8D22XB	Phadebas	Negative	Presumptive
8WXUW9	RSID	Inconclusive	Confirmatory
927828	Phadebas	Negative	Presumptive
AWM2QV	Phadebas	Negative	Presumptive
BLJ79Q	RSID	Negative	Confirmatory
BUXV7Y	Phadebas	Negative	Presumptive
CT4Y9C	Phadebas	Negative	Presumptive
CUWH34	Phadebas Press test	Negative	Presumptive
CWLNZ	RSID	Positive	Confirmatory
CY79E6	Phadebas	Negative	Presumptive
D8Y9K2	Phadebas	Inconclusive	Presumptive
DBGZH6	Phadebas	Negative	Presumptive
EBXDL8	Phadebas tube test	Negative	Presumptive
EXNXE3	Phadebas	Negative	Presumptive
F92HC4	Phadebas	Negative	Presumptive
FB83Y4	amylase diffusion	Negative	Presumptive
	micro (NEC)	Negative	Presumptive
FFJD6X	Phadebas	Inconclusive	Presumptive
FGBE6U	Phadebas	Negative	Presumptive
FPTPHV	Phadebas	Negative	Presumptive
GE84NK	RSID - Saliva	Negative	Confirmatory

TABLE 1 - Item 3

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
GZ76ZN	Rsid Saliva	Negative	Presumptive
H7PRBJ	RSID	Positive	Confirmatory
HAPXUW	Phadebas	Negative	Presumptive
HJCUHX	Phadebas	Negative	Presumptive
K28BTQ	Phadebas	Negative	Presumptive
K3WX4T	Phadebas	Inconclusive	Presumptive
KEA3BP	Phadebas	Negative	Presumptive
KGHKK2	α -Amylase	Negative	Confirmatory
KJN692	ALS	Positive	Presumptive
	RSID - Saliva	Negative	Confirmatory
L9N8AL	Phadebas	Negative	Presumptive
LQLPYG	RSID-Saliva	Negative	Confirmatory
M3B3FN	RSID - Saliva	Negative	Confirmatory
M48Q7U	Phadebas	Negative	Presumptive
MGAGMT	Phadebas	Negative	Presumptive
	Amilase[sic] - RSID	Negative	Confirmatory
N4Y7VW	Amylase Gel	Negative	Presumptive
	Micro	Negative	Presumptive
NALRPG	Phadebas Amylase	Negative	Presumptive
PF8YRL	Phadebas	Negative	Presumptive
R7EA6H	Phadebas	Negative	Presumptive
R7YFFN	Phadebas	Negative	Presumptive
TGBTYH	Salivary Amylase	Positive	Presumptive
ULKDJ9	ALS	Positive	Presumptive
	RSID-Saliva	Negative	Confirmatory
VK47JR	Phadebas	Negative	Presumptive
WAZNX9	RSID saliva	Negative	Presumptive
WLAQ8R	RSID Saliva	Negative	Confirmatory

TABLE 1 - Item 3

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
WQNYJB	Phadebas	Negative	Presumptive
Y7WPEE	Phadebas	Inconclusive	Presumptive
ZG8GHE	Phadebas	Negative	Presumptive
ZHYHGB	Phadebas	Negative	Presumptive
ZW33KZ	RSID-saliva	Positive	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
Phadebas	39	Negative
RSID-Saliva	14	Negative

TABLE 1 - Item 3

Other Body Fluid Types			
Webcode	Body Fluid Type - Test Performed	Results	Presumptive/Confirmatory
F92HC4	Semen - Choline/Florence Iodine	Positive	Presumptive
KGHKK2	Semen - ALS	Positive	Presumptive
MGAGMT	ALS	Positive	Presumptive
	Cells - Micro	Negative	Confirmatory
NALRPG	Urine - Creatinine	Negative	Presumptive
	Urine - DMAC	Negative	Presumptive
WLAQ8R	Urine - RSID Urine	Negative	Confirmatory

TABLE 1 - Item 4

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
2XTJD9	Otol	Positive	Presumptive
	Hemastix	Positive	Presumptive
4F8HWC	TMB	Positive	Presumptive
	Haematrace[sic]	Positive	Confirmatory
6DCQBC	KM	Positive	Presumptive
6RBUDH	KM	Positive	Presumptive
	Takayama	Positive	Confirmatory
	Ouchterlony-human	Positive	Confirmatory
7A8K9C	KM	Positive	Presumptive
7N9F6H	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7T4D8A	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7Y6RN4	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
8D22XB	KM	Positive	Presumptive
8WXUW9	visual	Positive	Presumptive
	Otol	Positive	Presumptive
927828	KM	Positive	Presumptive
AWM2QV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BLJ79Q	O-tol	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BUXV7Y	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
C74XV9	KM	Positive	Presumptive
CT4Y9C	KM	Positive	Presumptive

TABLE 1 - Item 4

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
CUWH34	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
CWLNMZ	O-tol	Positive	Presumptive
	RSID	Positive	Confirmatory
CY79E6	KM	Positive	Presumptive
D8Y9K2	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
DBGZH6	KM	Positive	Presumptive
DBYGLQ	KM	Positive	Presumptive
EBXDL8	KM	Positive	Presumptive
ELA9PQ	Visual	Positive	Presumptive
	LMG	Positive	Presumptive
	ABA HemaTrace	Positive	Confirmatory
EXNXE3	KM	Positive	Presumptive
F92HC4	kastle mayer[sic]	Positive	Presumptive
FB83Y4	o-tol	Positive	Presumptive
	hematrace	Positive	Presumptive
FFJD6X	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
FGBE6U	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
FPTPHV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
GE84NK	OT	Positive	Presumptive
	HT	Positive	Confirmatory
GUWFC9	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
GZ76ZN	Visual	Positive	Presumptive
	KM	Positive	Presumptive

TABLE 1 - Item 4

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
H4ZKQM	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
	Visual	Positive	Presumptive	
H7PRBJ	o-Tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
HAPXUW	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
HJCUHX	KM	Positive	Presumptive	
J8DD62	KM	Positive	Presumptive	
JVLQY7	KM	Positive	Presumptive	
K28BTQ	KM	Positive	Presumptive	
	ABACard Hematrace	Positive	Confirmatory	
K3WX4T	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
KEA3BP	KM	Positive	Presumptive	
	ABACard Hematrace	Positive	Confirmatory	
KETDZ6	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
KGHKK2	KM	Positive	Presumptive	
KJN692	Characteristic Appearance	Positive	Presumptive	
	KM	Positive	Presumptive	
L9N8AL	KM	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
LQLPYG	O-tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
LRV4H3	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
LUKB3Z	KM	Positive	Presumptive	

TABLE 1 - Item 4

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
M3B3FN	TMB (hemastix)	Positive	Presumptive
	Hematrace	Positive	Confirmatory
M48Q7U	KM	Positive	Presumptive
MGAGMT	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
N4Y7VW	O-Tol	Positive	Presumptive
	Hematrace	Positive	Presumptive
NALRPG	KM	Positive	Presumptive
	Micro (Mod Tak)	Positive	Confirmatory
PF8YRL	TMB	Positive	Presumptive
R7EA6H	KM	Positive	Presumptive
	ABAcad Hematrace	Positive	Confirmatory
R7YFFN	KM	Positive	Presumptive
RJZCCV	LMG	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
T9XLJE	Hemochromogen	Positive	Confirmatory
TGBTYH	O-tol	Positive	Presumptive
	Hematrace	Positive	Confirmatory
TGTPMB	KM	Positive	Presumptive
	Hematrace	Positive	Presumptive
ULKDJ9	o-Tol	Positive	Presumptive
	Hematrace	Positive	Confirmatory
VK47JR	KM	Positive	Presumptive
WAZNX9	Tetrabase test	Positive	Confirmatory
	RSID blood	Positive	Presumptive
WLAQ8R	Hemastix	Positive	Presumptive
	Hexagon OBTI	Positive	Confirmatory

TABLE 1 - Item 4

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
WQ4H49	Visual	Positive	Presumptive	
	LMG	Positive	Presumptive	
	ABAcad Hematrace	Positive	Confirmatory	
WQNYJB	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
XT8BX7	KM	Positive	Presumptive	
XXLHBR	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
Y7WPEE	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
ZG8GHE	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
ZHYHGB	Hematrace	Positive	Confirmatory	
	KM	Positive	Presumptive	
ZVKJNN	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
ZW33KZ	O-tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
KM	47	Positive
Hematrace	41	Positive

TABLE 1 - Item 4

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
2XTJD9	ALS - Polilight	Negative	Presumptive
	ALS - Crime lite	Negative	Presumptive
4F8HWC	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
6DCQBC	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
6RBUDH	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
7A8K9C	AP	Negative	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
7Y6RN4	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
8D22XB	AP	Negative	Presumptive
	P30 PSA	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
8WXUW9	ap	Negative	Presumptive
	micro	Negative	Confirmatory
927828	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
AWM2QV	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory

TABLE 1 - Item 4

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
BLJ79Q	AP	Inconclusive	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
BUXV7Y	AP	Negative	Presumptive	
	ABACard P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
C74XV9	AP	Negative	Presumptive	
CT4Y9C	AP	Negative	Presumptive	
	PSA	Negative	Presumptive	
	Micro (KPIC)	Negative	Confirmatory	
CWLNMZ	AP	Negative	Presumptive	
	RSID	Negative	Confirmatory	
CY79E6	AP	Negative	Presumptive	
	PSA (ABACard, P30)	Negative	Presumptive	
	Micro (KPIC)	Negative	Confirmatory	
D8Y9K2	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
DBGZH6	AP	Inconclusive	Presumptive	
	PSA (p30)	Negative	Presumptive	
	Microscopic	Negative	Confirmatory	
EBXDL8	Acid Phosphatase (AP)	Negative	Presumptive	
ELA9PQ	AP	Negative	Presumptive	
	ABA P30	Negative	Confirmatory	
	Microscope	Negative	Confirmatory	
EXNXE3	AP	Inconclusive	Presumptive	
	PSA	Negative	Presumptive	
	Micro	Negative	Confirmatory	
F92HC4	Acid Phosphatase	Negative	Presumptive	
	Choline/Florence Iodine	Negative	Presumptive	

TABLE 1 - Item 4

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
FB83Y4	Acid phosphatase	Negative	Presumptive	
	micro (sperm)	Negative	Confirmatory	
FFJD6X	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
FGBE6U	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
FPTPHV	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
GE84NK	ALS	Negative	Presumptive	
	AP	Negative	Presumptive	
GUWFC9	Micro	Negative	Confirmatory	
	PSA	Positive	Confirmatory	
GZ76ZN	AP	Negative		
	Micro	Negative		
	Rsid semen	Negative		
H4ZKQM	AP	Negative	Presumptive	
	PSA	Negative	Confirmatory	
	Micro	Negative	Confirmatory	
H7PRBJ	AP	Inconclusive	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
HAPXUW	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
HJCUHX	AP	Negative	Presumptive	
	P30/PSA	Negative	Presumptive	
	KPIC/Micro	Negative	Confirmatory	
J8DD62	AP	Negative	Presumptive	
JVLQY7	AP	Negative	Presumptive	

TABLE 1 - Item 4

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
K28BTQ	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
K3WX4T	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
KEA3BP	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
KETDZ6	Micro	Negative	Confirmatory
	P30	Positive	Confirmatory
KGHKK2	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
KJN692	ALS	Negative	Presumptive
	AP	Negative	Presumptive
L9N8AL	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
LQLPYG	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory
LRV4H3	Microscopic	Negative	Confirmatory
	PSA	Negative	Confirmatory
LUKB3Z	AP	Negative	Presumptive
M3B3FN	Micro	Negative	Confirmatory
	AP	Negative	Presumptive
M48Q7U	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro - KPIC	Negative	Confirmatory
MGAGMT	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
	Micro (CT)	Negative	Confirmatory

TABLE 1 - Item 4

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
N4Y7VW	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
NALRPG	AP	Negative	Presumptive
PF8YRL	Microscopy	Negative	Confirmatory
	AP	Negative	Presumptive
R7EA6H	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
	Micro	Negative	Confirmatory
R7YFFN	AP	Negative	Presumptive
	p30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
RJZCCV	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory
T9XLJE	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Positive	Presumptive
TGBTYH	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
ULKDJ9	ALS	Negative	Presumptive
	AP	Negative	Presumptive
	ABAcad P30	Negative	Confirmatory
	Sperm Cell Extraction/Micro	Negative	Confirmatory
VK47JR	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
WAZNX9	Microscopic analysis	Negative	Confirmatory
	PSA test	Negative	Presumptive
	RSID semen	Negative	Presumptive
WLAQ8R	ALS (Crime lite)	Negative	Presumptive
	Phosphatesmo	Negative	Confirmatory

TABLE 1 - Item 4

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
WQ4H49	Acid Phosphatase	Negative	Presumptive
	ABACard P30	Negative	Confirmatory
	Microscopic (sperm)	Negative	Confirmatory
WQNYJB	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
XXLHBR	Micro.	Negative	Confirmatory
	PSA	Negative	Confirmatory
Y7WPEE	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
ZG8GHE	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
ZHYHGB	ABACard P30	Negative	Presumptive
	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
ZW33KZ	AP	Inconclusive	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
AP	59	Negative
Microscopy	52	Negative

TABLE 1 - Item 4

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	Phadebas	Negative	Presumptive
4F8HWC	Phadebas amylase test	Negative	Presumptive
6DCQBC	Phadebas	Negative	Presumptive
6RBUDH	Phadebas	Negative	Presumptive
7A8K9C	Phadebas	Negative	Presumptive
7Y6RN4	Phadebas	Negative	Presumptive
8D22XB	Phadebas	Negative	Presumptive
8WXUW9	RSID	Negative	Confirmatory
927828	Phadebas	Negative	Presumptive
AWM2QV	Phadebas	Negative	Presumptive
BLJ79Q	RSID	Negative	Confirmatory
BUXV7Y	Phadebas	Negative	Presumptive
CT4Y9C	Phadebas	Negative	Presumptive
CWLMNZ	RSID	Positive	Confirmatory
CY79E6	Phadebas	Negative	Presumptive
D8Y9K2	Phadebas	Negative	Presumptive
DBGZH6	Phadebas	Negative	Presumptive
EBXDL8	Phadebas tube test	Negative	Presumptive
EXNXE3	Phadebas	Negative	Presumptive
F92HC4	Phadebas	Inconclusive	Presumptive
FB83Y4	amylase diffusion	Negative	Presumptive
	micro (NEC)	Negative	Presumptive
FFJD6X	Phadebas	Inconclusive	Presumptive
FGBE6U	Phadebas	Negative	Presumptive
FPTPHV	Phadebas	Negative	Presumptive
GE84NK	RSIS[sic] - Saliva	Negative	Confirmatory
GZ76ZN	rsid saliva	Negative	Presumptive
H7PRBJ	RSID	Negative	Confirmatory

TABLE 1 - Item 4

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
HAPXUW	Phadebas	Negative	Presumptive
HJCUHX	Phadebas	Negative	Presumptive
K28BTQ	Phadebas	Negative	Presumptive
K3WX4T	Phadebas	Negative	Presumptive
KEA3BP	Phadebas	Negative	Presumptive
KGHKK2	α-Amylase	Negative	Confirmatory
KJN692	ALS	Negative	Presumptive
	RSID - Saliva	Negative	Confirmatory
L9N8AL	Phadebas	Negative	Presumptive
LQLPYG	RSID-Saliva	Negative	Confirmatory
M3B3FN	RSID - Saliva	Negative	Confirmatory
M48Q7U	Phadebas	Negative	Presumptive
MGAGMT	Phadebas	Negative	Presumptive
	Amilase[sic] - RSID	Negative	Confirmatory
N4Y7WW	Amylase Gel	Negative	Presumptive
	Micro	Negative	Presumptive
NALRPG	Phadebas Amylase	Negative	Presumptive
PF8YRL	Phadebas	Negative	Presumptive
R7EA6H	Phadebas	Negative	Presumptive
R7YFFN	Phadebas	Negative	Presumptive
TGBTYH	Salivary amylase	Negative	Presumptive
ULKDJ9	ALS	Negative	Presumptive
	RSID-Saliva	Negative	Confirmatory
VK47JR	Phadebas	Negative	Presumptive
WAZNX9	RSID Saliva	Negative	Presumptive
WLAQ8R	RSID Saliva	Negative	Confirmatory
WQNYJB	Phadebas	Negative	Presumptive
Y7WPEE	Phadebas	Inconclusive	Presumptive

TABLE 1 - Item 4

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
ZG8GHE	Phadebas	Negative	Presumptive
ZHYHGB	Phadebas	Negative	Presumptive
ZW33KZ	RSID-Saliva	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
Phadebas	37	Negative
RSID-Saliva	13	Negative

TABLE 1 - Item 4

Other Body Fluid Types			
Webcode	Body Fluid Type - Test Performed	Results	Presumptive/Confirmatory
KGHKK2	Semen - ALS	Negative	Presumptive
MGAGMT	ALS	Negative	
	Cells - Micro	Negative	

TABLE 1 - Item 5

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
2XTJD9	Otol	Positive	Presumptive
	Hemastix	Positive	Presumptive
4F8HWC	TMB	Positive	Presumptive
	Haematrace[sic]	Positive	Confirmatory
6DCQBC	KM	Positive	Presumptive
6RBUDH	KM	Positive	Presumptive
	Takayama	Positive	Confirmatory
	Ouchterlony-human	Positive	Confirmatory
7A8K9C	KM	Positive	Presumptive
7N9F6H	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7T4D8A	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7Y6RN4	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
8D22XB	KM	Positive	Presumptive
8WXUW9	visual	Positive	Presumptive
	Otol	Positive	Presumptive
927828	KM	Positive	Presumptive
AWM2QV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BLJ79Q	O-tol	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BUXV7Y	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
C74XV9	KM	Positive	Presumptive
CT4Y9C	KM	Positive	Presumptive

TABLE 1 - Item 5

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
CUWH34	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
CWLNMZ	O-tol	Positive	Presumptive
	RSID	Positive	Confirmatory
CY79E6	KM	Positive	Presumptive
D8Y9K2	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
DBGZH6	KM	Positive	Presumptive
DBYGLQ	KM	Positive	Presumptive
EBXDL8	KM	Positive	Presumptive
ELA9PQ	Visual	Positive	Presumptive
	LMG	Positive	Presumptive
	ABA HemaTrace	Positive	Confirmatory
EXNXE3	KM	Positive	Presumptive
F92HC4	kastle mayer[sic]	Positive	Presumptive
FB83Y4	o-tol	Positive	Presumptive
	hematrace	Positive	Presumptive
	IR light & camera	Negative	Presumptive
FFJD6X	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
FGBE6U	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
FPTPHV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
GE84NK	OT	Positive	Presumptive
	HT	Positive	Confirmatory
GUWFC9	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory

TABLE 1 - Item 5

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
GZ76ZN	Visual, ALS	Positive	Presumptive	
	KM	Positive	Presumptive	
	Haematrace[sic]	Positive	Presumptive	
	low power micro	Positive	Presumptive	
H4ZKQM	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
	Visual	Positive	Presumptive	
H7PRBJ	o-Tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
HAPXUW	KM	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
HJCUHX	KM	Positive	Presumptive	
J8DD62	KM	Positive	Presumptive	
JVLQY7	KM	Positive	Presumptive	
K28BTQ	KM	Positive	Presumptive	
	ABAcad Hematrace	Positive	Confirmatory	
K3WX4T	KM	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
KEA3BP	KM	Positive	Presumptive	
	ABAcad Hematrace	Positive	Confirmatory	
KETDZ6	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
KGHKK2	KM	Positive	Presumptive	
KJN692	Characteristic Appearance	Positive	Presumptive	
	KM	Positive	Presumptive	
L9N8AL	KM	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
LQLPYG	O-tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	

TABLE 1 - Item 5

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
LRV4H3	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
LUKB3Z	KM	Positive	Presumptive	
M3B3FN	TMB (hemastix)	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
M48Q7U	KM	Positive	Presumptive	
MGAGMT	KM	Positive	Presumptive	
	Henatrace[sic]	Positive	Confirmatory	
N4Y7VW	O-Tol	Positive	Presumptive	
	Hematrace	Positive	Presumptive	
NALRPG	KM	Positive	Presumptive	
	Micro (Mod Tak)	Positive	Confirmatory	
PF8YRL	TMB	Positive	Presumptive	
R7EA6H	KM	Positive	Presumptive	
	ABACard Hematrace	Positive	Confirmatory	
R7YFFN	KM	Positive	Presumptive	
RJZCCV	LMG	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
T9XLJE	Hemastix (TMB)	Positive	Presumptive	
	Hemochromogen	Positive	Confirmatory	
TGPTYH	O-tol	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
TGTPMB	KM	Positive	Presumptive	
	Hematrace	Positive	Presumptive	
ULKDJ9	o-Tol	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
VK47JR	KM	Positive	Presumptive	

TABLE 1 - Item 5

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
WAZNX9	Oblique light + binoculair[sic]	Positive	Presumptive
	Tetrabase test	Positive	Confirmatory
	RSID blood	Positive	Presumptive
WLAQ8R	Hemastix	Positive	Presumptive
	Hexagon OBTI	Positive	Confirmatory
WQ4H49	Visual	Positive	Presumptive
	LMG	Positive	Presumptive
	ABAcad Hematrace	Positive	Confirmatory
WQNYJB	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
XT8BX7	KM	Positive	Presumptive
Y7WPEE	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
ZG8GHE	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
ZHYHGB	Hematrace	Positive	Confirmatory
	KM	Positive	Presumptive
ZVKJNN	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
ZW33KZ	O-tol	Positive	Presumptive
	HemaTrace	Positive	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
KM	47	Positive
Hematrace	40	Positive

TABLE 1 - Item 5

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	AP	Negative	Presumptive
	ABACard P30	Negative	Presumptive
2XTJD9	ALS - Polilight	Negative	Presumptive
	ALS - Crime lite	Negative	Presumptive
4F8HWC	AP	Negative	Presumptive
6DCQBC	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
6RBUDH	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	ABACard	Negative	Confirmatory
7A8K9C	AP	Negative	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
7N9F6H	ALS	Negative	Presumptive
7Y6RN4	AP	Negative	Presumptive
	ABACard P30	Negative	Presumptive
	Micro		Confirmatory
8D22XB	AP	Negative	Presumptive
	P30 PSA	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
8WXUW9	N/A given scenario		
927828	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
AWM2QV	AP	Negative	Presumptive
	ABACard P30	Negative	Presumptive
BLJ79Q	ALS	Negative	Presumptive
	AP	Negative	Presumptive

TABLE 1 - Item 5

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
BUXV7Y	AP	Negative	Presumptive	
	ABACard P30	Negative	Presumptive	
CT4Y9C	AP	Negative	Presumptive	
	PSA	Negative	Presumptive	
	Micro (KPIC)	Negative	Confirmatory	
CWLNMZ	AP	Negative	Presumptive	
	RSID	Negative	Confirmatory	
CY79E6	AP	Negative	Presumptive	
	PSA (ABACard, P30)	Negative	Presumptive	
	Micro (KPIC)	Negative	Confirmatory	
DBGZH6	AP	Inconclusive	Presumptive	
	PSA (p30)	Negative	Presumptive	
	Microscopic	Negative	Confirmatory	
	ALS	Negative	Presumptive	
DBYGLQ	AP	Negative		
EBXDL8	Acid Phosphatase (AP) (Blotting method)	Negative	Presumptive	
ELA9PQ	ALS	Positive	Presumptive	
	AP	Negative	Presumptive	
	ABA P30	Negative	Confirmatory	
	Microscope	Negative	Confirmatory	
EXNXE3	AP	Negative	Presumptive	
	PSA	Negative	Presumptive	
	Micro	Negative	Confirmatory	
FB83Y4	Acid phosphatase	Negative	Presumptive	
	micro (sperm)	Negative	Confirmatory	
	ALS	Negative	Presumptive	
FGBE6U	AP	Negative	Presumptive	
FPTPHV	AP	Negative	Presumptive	
	ABACard P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	

TABLE 1 - Item 5

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
GE84NK	ALS	Negative	Presumptive
	AP	Negative	Presumptive
GZ76ZN	AP	Negative	Presumptive
	micro on blood stain	Negative	Confirmatory
H4ZKQM	ALS	Negative	Presumptive
	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
	Micro	Negative	Confirmatory
H7PRBJ	ALS	Positive	Presumptive
	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory
HJCUHX	AP	Negative	Presumptive
	P30/PSA	Negative	Presumptive
	KPIC/Micro	Negative	Confirmatory
	ALS	Negative	Presumptive
K28BTQ	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
KEA3BP	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
KGHKK2	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
KJN692	ALS	Negative	Presumptive
	AP	Negative	Presumptive
L9N8AL	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
LQLPYG	AP	Negative	Presumptive
	ALS	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory

TABLE 1 - Item 5

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
M3B3FN	AP	Negative	Presumptive
M48Q7U	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro - KPIC	Negative	Confirmatory
MGAGMT	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
	Micro (CT)	Negative	Confirmatory
N4Y7VW	ALS	Negative	Presumptive
	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Presumptive
NALRPG	ALS	Negative	
PF8YRL	AP	Negative	Presumptive
R7EA6H	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
R7YFFN	ALS	Negative	Presumptive
	AP	Negative	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
RJZCCV	ALS	Negative	Presumptive
TGBTYH	AP	Negative	Presumptive
ULKDJ9	ALS	Negative	Presumptive
	AP	Negative	Presumptive
	ABAcadP30	Negative	Confirmatory
	Sperm Cell Extraction/Micro	Negative	Confirmatory
VK47JR	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory

TABLE 1 - Item 5

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
WAZNX9	Crimescope (css/455nm)	Negative	Presumptive
	PSA test	Negative	Presumptive
	RSID semen	Negative	Presumptive
	Microscopic analysis	Negative	
WLAQ8R	ALS (Crime lite)	Negative	Presumptive
WQ4H49	ALS	Positive	Presumptive
	Acid Phosphatase	Negative	Presumptive
	ABAcad P30	Negative	Confirmatory
	Microscopic (sperm)	Negative	Confirmatory
ZHYHGB	ABAcad P30	Negative	Presumptive
	AP	Negative	Presumptive
ZVKJNN	ALS	Negative	Presumptive
ZW33KZ	ALS	Negative	Presumptive
	AP	Negative	Presumptive
	micro	Negative	Confirmatory
	PSA	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
AP	46	Negative
P30/PSA	33	Negative

TABLE 1 - Item 5

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	Phadebas	Negative	Presumptive
4F8HWC	Phadebas amylase test (paper)	Negative	Presumptive
6DCQBC	Phadebas	Negative	Presumptive
7A8K9C	Phadebas	Negative	Presumptive
7T4D8A	Phadebas Press Test	Negative	Presumptive
	Phadebas Tube Test on Blood	Inconclusive	Presumptive
7Y6RN4	Phadebas	Negative	Presumptive
8D22XB	Phadebas	Negative	Presumptive
8WXUW9	ALS	Negative	Presumptive
927828	Phadebas	Negative	Presumptive
AWM2QV	Phadebas	Negative	Presumptive
BUXV7Y	Phadebas	Negative	Presumptive
CT4Y9C	Phadebas	Negative	Presumptive
CUWH34	Phadebas press test	Negative	Presumptive
	Phadebas tube test	Inconclusive	Presumptive
CWLNMZ	RSID	Positive	Confirmatory
CY79E6	Phadebas	Negative	Presumptive
D8Y9K2	Phadebas	Negative	Presumptive
DBGZH6	Phadebas	Negative	Presumptive
	ALS	Negative	Presumptive
EBXDL8	Phadebas Paper (blotting method)	Negative	Presumptive
EXNXE3	Phadebas	Negative	Presumptive
FB83Y4	amylase diffusion	Negative	Presumptive
	micro (NEC)	Negative	Presumptive
FFJD6X	Phadebas	Negative	Presumptive
FGBE6U	Phadebas	Negative	Presumptive
FPTPHV	Phadebas	Negative	Presumptive
GE84NK	RSID - Saliva	Negative	Confirmatory

TABLE 1 - Item 5

Saliva				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
GZ76ZN	Phadebas paper test	Negative	Presumptive	
H7PRBJ	RSID	Negative	Confirmatory	
HAPXUW	Phadebas	Negative	Presumptive	
HJCUHX	Phadebas	Negative	Presumptive	
K28BTQ	Phadebas	Negative	Presumptive	
K3WX4T	Phadebas	Negative	Presumptive	
KEA3BP	Phadebas	Negative	Presumptive	
KGHKK2	α Amylase	Negative	Confirmatory	
KJN692	ALS	Negative	Presumptive	
	RSID Saliva	Negative	Confirmatory	
L9N8AL	Phadebas	Negative	Presumptive	
LQLPYG	RSID-Saliva	Negative	Confirmatory	
M48Q7U	Phadebas	Negative	Presumptive	
MGAGMT	Phadebas	Negative	Presumptive	
	Amilase[sic] - RSID	Negative	Confirmatory	
N4Y7VW	ALS	Negative	Presumptive	
	Amylase Gel	Negative	Presumptive	
	Micro	Negative	Presumptive	
NALRPG	ALS	Negative		
PF8YRL	Phadebas	Negative	Presumptive	
R7EA6H	Phadebas	Negative	Presumptive	
R7YFFN	Phadebas	Negative	Presumptive	
ULKDJ9	ALS	Negative	Presumptive	
	RSID-Saliva	Negative	Confirmatory	
VK47JR	Phadebas	Negative	Presumptive	
WAZNX9	Crimescope (css/455 nm)	Negative	Presumptive	
	RSID saliva	Negative	Presumptive	
WLAQ8R	RSID Saliva	Negative	Confirmatory	

TABLE 1 - Item 5

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
WQNYJB	Phadebas	Negative	Presumptive
Y7WPEE	Phadebas	Negative	Presumptive
ZG8GHE	Phadebas	Negative	Presumptive
ZHYHGB	Phadebas	Negative	Presumptive
ZVKJNN	ALS	Negative	Presumptive
ZW33KZ	RSID-saliva	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
Phadebas	37	Negative
RSID-Saliva	10	Negative

TABLE 1 - Item 5

Other Body Fluid Types			
Webcode	Body Fluid Type - Test Performed	Results	Presumptive/Confirmatory
6DCQBC	Semen/Saliva - ALS	Negative	Presumptive
7A8K9C	Semen/Saliva - ALS	Negative	Presumptive
7N9F6H	epithelial cells - swabbed for possible e.c.'s		Presumptive
8D22XB	Semen/Saliva - ALS	Negative	Presumptive
927828	ALS	Negative	Presumptive
CT4Y9C	Semen/Saliva - ALS	Negative	Presumptive
CWLNMZ	Semen - Polilight	Negative	Presumptive
	Saliva - Polilight	Negative	Presumptive
CY79E6	Semen, Saliva - ALS	Negative	Presumptive
D8Y9K2	Semen - ALS	Negative	Presumptive
DBYGLQ	ALS		
GZ76ZN	Epithelial - the inside and the outside of glove was minitape separately for epithelial cells		
KGHKK2	Semen - ALS	Negative	Presumptive
LRV4H3	Epithelial cells - Microscopic slide not made to conserve sample. Glove was swabbed.		
M3B3FN	Saliva/Semen - ALS (polilight)	Negative	Presumptive
M48Q7U	Semen/Saliva - ALS	Negative	Presumptive
MGAGMT	ALS	Negative	Presumptive
T9XLJE	ALS	Negative	Presumptive
VK47JR	Semen - ALS	Negative	Presumptive
WQNYJB	Semen - ALS	Negative	Presumptive
XT8BX7	Semen/Saliva - ALS (Alternate light source)	Negative	Presumptive
XXLHBR	Epithelial cells - ALS *see comments section		

TABLE 1 - Item 6

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
2XTJD9	Otol	Positive	Presumptive
	Hemastix	Positive	Presumptive
4F8HWC	TMB	Positive	Presumptive
	Haematrace[sic]	Positive	Confirmatory
6DCQBC	KM	Positive	Presumptive
6RBUDH	KM	Positive	Presumptive
	Takayama	Positive	Confirmatory
	Ouchterlony-human	Positive	Confirmatory
7A8K9C	KM	Positive	Presumptive
7N9F6H	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7T4D8A	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
7Y6RN4	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
8D22XB	KM	Positive	Presumptive
8WXUW9	visual	Positive	Presumptive
	otol	Positive	Presumptive
927828	KM	Positive	Presumptive
AWM2QV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BLJ79Q	O-tol	Positive	Presumptive
	Hematrace	Positive	Confirmatory
BUXV7Y	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
C74XV9	KM	Positive	Presumptive
CT4Y9C	KM	Positive	Presumptive

TABLE 1 - Item 6

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
CUWH34	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
CWLNMZ	O-tol	Positive	Presumptive
	RSID	Positive	Confirmatory
CY79E6	KM	Positive	Presumptive
D8Y9K2	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
DBGZH6	KM	Positive	Presumptive
DBYGLQ	KM	Positive	Presumptive
EBXDL8	KM	Positive	Presumptive
ELA9PQ	visual	Positive	Presumptive
	LMG	Positive	Presumptive
	ABA Hema Trace	Positive	Confirmatory
EXNXE3	KM (Q1)	Negative	Presumptive
	KM (Q2)	Negative	Presumptive
	KM (Q3)	Positive	Presumptive
F92HC4	kastle mayer[sic]	Positive	Presumptive
	Hematrace	Positive	Confirmatory
FB83Y4	o-tol	Positive	Presumptive
	hematrace	Positive	Presumptive
FFJD6X	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
FGBE6U	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
FPTPHV	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory
GE84NK	OT	Positive	Presumptive
	HT	Positive	Confirmatory

TABLE 1 - Item 6

Blood			
Webcode	Test Performed	Results	Presumptive/Confirmatory
GUWFC9	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
GZ76ZN	Visual	Positive	Presumptive
	Km	Positive	Presumptive
	Haematrace[sic]	Positive	Presumptive
H4ZKQM	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
	Visual	Positive	Presumptive
H7PRBJ	o-Tol	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
HAPXUW	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
HJCUHX	KM	Positive	Presumptive
J8DD62	KM	Positive	Presumptive
JVLQY7	KM	Positive	Presumptive
K28BTQ	KM	Positive	Presumptive
	ABAcad Hematrace	Positive	Confirmatory
K3WX4T	KM	Positive	Presumptive
	HemaTrace	Positive	Confirmatory
KEA3BP	KM	Positive	Presumptive
	ABAcad Hematrace	Positive	Confirmatory
KETDZ6	LMG	Positive	Presumptive
	Hematrace	Positive	Confirmatory
KGHKK2	KM	Positive	Presumptive
KJN692	Characteristic Appearance	Positive	Presumptive
	KM	Positive	Presumptive
L9N8AL	KM	Positive	Presumptive
	Hematrace	Positive	Confirmatory

TABLE 1 - Item 6

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
LQLPYG	O-tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
LRV4H3	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
LUKB3Z	KM	Positive	Presumptive	
M3B3FN	TMB (hemastix)	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
M48Q7U	KM	Positive	Presumptive	
MGAGMT	KM	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
N4Y7VW	O-Tol	Positive	Presumptive	
	Hematrace	Positive	Presumptive	
NALRPG	KM	Positive	Presumptive	
	Micro (Mod Tak)	Positive	Confirmatory	
PF8YRL	TMB	Positive	Presumptive	
R7EA6H	KM	Positive	Presumptive	
	ABACard Hematrace	Positive	Confirmatory	
R7YFFN	KM	Positive	Presumptive	
RJZCCV	LMG	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
T9XLJE	Hemochromogen	Positive	Confirmatory	
TGBTYH	O-tol	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
TGTPMB	KM	Positive	Presumptive	
	Hematrace	Positive	Presumptive	
ULKDJ9	o-Tol	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
VK47JR	KM	Positive	Presumptive	

TABLE 1 - Item 6

Blood				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
WAZNX9	Oblique light	Positive	Presumptive	
	Tetrabase test	Positive	Confirmatory	
	RSID blood	Positive	Presumptive	
WLAQ8R	Hemastix	Positive	Presumptive	
	Hexagon OBTI	Positive	Confirmatory	
WQ4H49	Visual	Positive	Presumptive	
	LMG	Positive	Presumptive	
	ABAcad Hematrace	Positive	Confirmatory	
WQNYJB	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
XT8BX7	KM	Positive	Presumptive	
XXLHBR	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
Y7WPEE	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
ZG8GHE	KM	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	
ZHYHGB	Hematrace	Positive	Confirmatory	
	KM	Positive	Presumptive	
ZVKJNN	LMG	Positive	Presumptive	
	Hematrace	Positive	Confirmatory	
ZW33KZ	O-tol	Positive	Presumptive	
	HemaTrace	Positive	Confirmatory	

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
KM	47	Positive
Hematrace	45	Positive

TABLE 1 - Item 6

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	AP	Negative	Presumptive
	ABACard P30	Negative	Presumptive
2XTJD9	ALS - Polilight	Positive	Presumptive
	ALS - Crime lite	Positive	Presumptive
4F8HWC	AP	Negative	Presumptive
6DCQBC	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
6RBUDH	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
	ABACard	Negative	Confirmatory
7A8K9C	AP	Negative	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
7N9F6H	ALS	Negative	Presumptive
7Y6RN4	AP	Negative	Presumptive
	ABACard P30	Negative	Presumptive
	Micro	Negative	Confirmatory
8D22XB	AP	Negative	Presumptive
	P30 PSA	Negative	Presumptive
	KPIC (Micro)	Negative	Confirmatory
8WXUW9	ap	Negative	Presumptive
927828	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
AWM2QV	AP	Negative	Presumptive
	ABACard P30	Negative	Presumptive
BLJ79Q	ALS	Negative	Presumptive
	AP	Negative	Presumptive

TABLE 1 - Item 6

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
BUXV7Y	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
C74XV9	ALS	Negative	Presumptive
CT4Y9C	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro (KPIC)	Negative	Confirmatory
CUWH34	AP	Negative	Presumptive
	Micro	Negative	Confirmatory
CWLNMZ	AP	Negative	Presumptive
	RSID	Negative	Confirmatory
CY79E6	AP	Negative	Presumptive
	PSA (ABAcad, P30)	Negative	Presumptive
	Micro (KPIC)	Negative	Confirmatory
DBGZH6	AP	Inconclusive	Presumptive
	PSA (p30)	Negative	Presumptive
	Microscopic	Negative	Confirmatory
	ALS	Negative	Presumptive
DBYGLQ	AP	Negative	
EBXDL8	Acid Phosphatase (AP) (blotting method)	Negative	Presumptive
ELA9PQ	ALS	Positive	Presumptive
	AP	Negative	Presumptive
	ABA P30	Negative	Confirmatory
	Microscope	Negative	Confirmatory
EXNXE3	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
F92HC4	Acid phosphatase	Negative	Presumptive

TABLE 1 - Item 6

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
FB83Y4	Acid phosphatase	Negative	Presumptive	
	micro (sperm)	Negative	Confirmatory	
	ALS	Negative	Presumptive	
FGBE6U	AP	Negative	Presumptive	
FPTPHV	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
	Micro	Negative	Confirmatory	
GE84NK	ALS	Negative	Presumptive	
	AP	Negative	Presumptive	
GZ76ZN	i/s & o/s ap tested	Negative	Presumptive	
	micro on blood stain	Negative	Confirmatory	
H4ZKQM	ALS	Positive	Presumptive	
	AP	Negative	Presumptive	
	PSA	Negative		
	Micro	Negative		
H7PRBJ	ALS	Positive	Presumptive	
	AP	Negative	Presumptive	
	Micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	
HJCUHX	AP	Negative	Presumptive	
	P30/PSA	Negative	Presumptive	
	KPIC/Micro	Negative	Confirmatory	
	ALS	Negative	Presumptive	
J8DD62	AP	Negative	Presumptive	
	ALS	Negative	Presumptive	
K28BTQ	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	
KEA3BP	AP	Negative	Presumptive	
	ABAcad P30	Negative	Presumptive	

TABLE 1 - Item 6

Semen			
Webcode	Test Performed	Results	Presumptive/Confirmatory
KGHKK2	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
KJN692	ALS	Negative	Presumptive
	AP	Negative	Presumptive
L9N8AL	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
LQLPYG	AP	Negative	Presumptive
	ALS	Negative	Presumptive
	Micro	Negative	Confirmatory
	PSA	Negative	Confirmatory
LUKB3Z	ALS	Negative	Presumptive
M3B3FN	AP	Negative	Presumptive
M48Q7U	AP	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro - KPIC	Negative	Confirmatory
MGAGMT	AP	Negative	Presumptive
	PSA	Negative	Confirmatory
	Micro (CT)	Negative	Confirmatory
N4Y7VW	ALS	Negative	Presumptive
	PSA	Negative	Presumptive
	Micro	Negative	Confirmatory
NALRPG	ALS	Negative	
PF8YRL	AP	Negative	Presumptive
R7EA6H	AP	Negative	Presumptive
	ABAcad P30	Negative	Presumptive
R7YFFN	ALS	Negative	Presumptive
	AP	Negative	Presumptive
	P30 (PSA)	Negative	Presumptive
	KPIC (micro)	Negative	Confirmatory
RJZCCV	ALS	Negative	Presumptive

TABLE 1 - Item 6

Semen				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
TGBTYH	AP	Negative	Presumptive	
ULKDJ9	ALS	Negative	Presumptive	
	AP	Negative	Presumptive	
	ABAcad P30	Negative	Confirmatory	
	Sperm Cell Extraction/Micro	Negative	Confirmatory	
VK47JR	AP	Negative	Presumptive	
	PSA	Negative	Presumptive	
	Micro	Negative	Confirmatory	
WAZNX9	Crimescope (css/455 nm)	Negative	Presumptive	
	PSA test	Negative	Presumptive	
	RSID semen	Negative	Presumptive	
	Phosphatesmo test	Negative	Presumptive	
	Microscopic analysis	Negative		
WLAQ8R	ALS (Crime lite)	Negative	Presumptive	
WQ4H49	ALS	Negative	Presumptive	
	Acid Phosphatase	Negative	Presumptive	
	ABAcad P30	Negative	Confirmatory	
	Microscope (sperm)	Negative	Confirmatory	
ZHYHGB	ABAcad P30	Negative	Presumptive	
	AP	Negative	Presumptive	
ZVKJNN	ALS	Negative	Presumptive	
ZW33KZ	AP	Negative	Presumptive	
	ALS	Negative	Presumptive	
	micro	Negative	Confirmatory	
	PSA	Negative	Confirmatory	

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
AP	49	Negative
P30/PSA	33	Negative

TABLE 1 - Item 6

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
2NMLP6	Phadebas	Negative	Presumptive
4F8HWC	Phadebas amylase test (paper)	Negative	Presumptive
6DCQBC	Phadebas	Negative	Presumptive
7A8K9C	Phadebas	Negative	Presumptive
7T4D8A	Phadebas Press Test	Negative	Presumptive
	Phadebas Tube Test on Blood	Inconclusive	Presumptive
7Y6RN4	Phadebas	Negative	Presumptive
8D22XB	Phadebas	Negative	Presumptive
8WXUW9	RSID	Negative	Confirmatory
927828	Phadebas	Negative	Presumptive
AWM2QV	Phadebas	Negative	Presumptive
BUXV7Y	Phadebas	Negative	Presumptive
CT4Y9C	Phadebas	Negative	Presumptive
CUWH34	Phadebas tube test	Inconclusive	Presumptive
CWLNMZ	RSID	Positive	Confirmatory
CY79E6	Phadebas	Negative	Presumptive
D8Y9K2	Phadebas	Negative	Presumptive
DBGZH6	Phadebas	Negative	Presumptive
	ALS	Negative	Presumptive
EBXDL8	Phadebas Paper (blotting method)	Negative	Presumptive
EXNXE3	Phadebas	Negative	Presumptive
FB83Y4	amylase diffusion	Negative	Presumptive
	micro (NEC)	Negative	Presumptive
FFJD6X	Phadebas	Negative	Presumptive
FGBE6U	Phadebas	Negative	Presumptive
FPTPHV	Phadebas	Negative	Presumptive
GE84NK	RSID - Saliva	Negative	Confirmatory
GZ76ZN	phadebas paper test on o/s ft and back	Negative	Presumptive

TABLE 1 - Item 6

Saliva				
Webcode	Test Performed	Results	Presumptive/Confirmatory	
H7PRBJ	RSID	Negative	Confirmatory	
HAPXUW	Phadebas	Negative	Presumptive	
HJCUHX	Phadebas	Negative	Presumptive	
K28BTQ	Phadebas	Negative	Presumptive	
K3WX4T	Phadebas	Negative	Presumptive	
KEA3BP	Phadebas	Negative	Presumptive	
KGHKK2	α -Amylase	Negative	Confirmatory	
KJN692	ALS	Negative	Presumptive	
	RSID - Saliva	Negative	Confirmatory	
L9N8AL	Phadebas	Negative	Presumptive	
LQLPYG	RSID-Saliva	Negative	Confirmatory	
M48Q7U	Phadebas	Negative	Presumptive	
MGAGMT	Phadebas	Negative	Presumptive	
	Amilase[sic] - RSID	Negative	Confirmatory	
N4Y7WW	ALS	Negative	Presumptive	
	Amylase Gel	Negative	Presumptive	
	Micro	Negative	Presumptive	
NALRPG	ALS	Negative		
PF8YRL	Phadebas	Negative	Presumptive	
R7EA6H	Phadebas	Negative	Presumptive	
R7YFFN	Phadebas	Negative	Presumptive	
ULKDJ9	ALS	Negative	Presumptive	
	RSID-Saliva	Negative	Confirmatory	
VK47JR	Phadebas	Negative	Presumptive	
WAZNX9	Crimescope (css/455 nm)	Negative	Presumptive	
	RSID saliva	Negative	Presumptive	
WLAQ8R	RSID Saliva	Negative	Confirmatory	
WQNYJB	Phadebas	Negative	Presumptive	

TABLE 1 - Item 6

Saliva			
Webcode	Test Performed	Results	Presumptive/Confirmatory
Y7WPEE	Phadebas	Negative	Presumptive
ZG8GHE	Phadebas	Negative	Presumptive
ZHYHGB	Phadebas	Negative	Presumptive
ZVKJNN	ALS	Negative	Presumptive
ZW33KZ	RSID-saliva	Negative	Confirmatory

Response Summary		
<u>Test Performed</u>	<u>Total Reported</u>	<u>Consensus Result</u>
Phadebas	37	Negative
RSID-Saliva	11	Negative

TABLE 1 - Item 6

Other Body Fluid Types			
Webcode	Body Fluid Type - Test Performed	Results	Presumptive/Confirmatory
6DCQBC	Semen/Saliva - ALS	Negative	Presumptive
7A8K9C	Semen/Saliva - ALS	Negative	Presumptive
7N9F6H	epithelial cells - swabbed swabbed[sic] for poss. e.c.'s		Presumptive
8D22XB	Semen/Saliva - ALS	Negative	Presumptive
927828	ALS	Negative	Presumptive
CT4Y9C	Semen/Saliva - ALS	Negative	Presumptive
CWLNMZ	Semen - Polilight	Negative	Presumptive
	Saliva - Polilight	Negative	Presumptive
CY79E6	Semen, Saliva - ALS	Negative	Presumptive
D8Y9K2	Semen - ALS	Negative	Presumptive
DBYGLQ	ALS		
GZ76ZN	Fibres - Fibre lifted o/s ft and back of top		
KGHKK2	Semen - ALS	Negative	Presumptive
LRV4H3	Epithelial cells - Microscopic slide not made to conserve sample. Tank top was swabbed.		
M3B3FN	Saliva/Semen - ALS (Polilight)	Negative	Presumptive
M48Q7U	Semen/Saliva - ALS	Negative	Presumptive
MGAGMT	ALS	Negative	Presumptive
T9XLJE	ALS	Negative	Presumptive
VK47JR	Semen - ALS	Negative	Presumptive
WQNYJB	Semen - ALS	Negative	Presumptive
XT8BX7	Semen/Saliva - ALS	Negative	Presumptive
XXLHBR	Epithelial cells - *see comments section		

Screening Conclusions

Indicate the overall conclusion from the results of the screening tests performed.

TABLE 2 - Item 1
Four swabs spotted with 100µl of Saliva.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
2NMLP6	Negative	Negative	Positive
2XTJD9	Negative	Negative	
4F8HWC	Blood not detected	Semen not detected	Saliva detected*
6DCQBC	No blood indicated.	No semen identified.	Saliva indicated.
6RBUDH	Negative	Negative	Amylase positive
7A8K9C	No blood indicated.	No semen identified.	Saliva indicated.
7N9F6H	+ positive	- negative	possible
7T4D8A	Blood not detected	Semen not detected	Amylase indicative of saliva
7Y6RN4	Negative	Negative	Positive
8D22XB	No blood indicated	No semen identified	Saliva indicated
8WXUW9	Negative	Negative	Positive: saliva present
927828	No blood indicated	No semen identified	Saliva indicated
AWM2QV	Negative	Negative	Positive
BLJ79Q	Negative	Negative	Positive
BUXV7Y	Negative	Negative	Positive
C74XV9	Negative	Negative	Positive
CT4Y9C	No blood indicated.	No semen identified.	Saliva indicated.
CUWH34	Blood not detected	Semen not detected	Amylase, indicative of saliva detected
CWLNMZ	Negative	Positive	Positive

TABLE 2 - Item 1
Four swabs spotted with 100µl of Saliva.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
CWLNMZ	Negative	Positive	Positive
CY79E6	No blood indicated.	No semen identified.	Saliva indicated.
D8Y9K2	Blood not detected	Semen not detected	Amylase, indicative of saliva
DBGZH6	No blood indicated.	No semen identified.	Saliva indicated.
DBYGLQ	Negative	Negative	Positive
EBXDL8	Negative	Negative	Positive for amylase, a constituent of saliva
ELA9PQ	no blood indicated	no seminal fluid indicated	
EXNXE3	No blood indicated	No semen identified	Saliva indicated
F92HC4	No blood detected	No semen detected	Possible saliva detected
FB83Y4	Negative	Negative	Positive
FFJD6X	Blood not detected	Semen not detected	Amylase, indicative of saliva
FGBE6U	Negative	Negative	Positive
FPTPHV	Negative	Negative	Positive
GE84NK	Negative	Negative	Positive (see comments section)
GUWFC9	Positive	Negative	Possible
GZ76ZN	Negative	Negative	Positive: indicative for human saliva
H4ZKQM	no blood was indicated	no seminal fluid indicated	
H7PRBJ	Negative	Negative	Positive
HAPXUW	Blood not detected	Semen not detected	Amylase, indicative of saliva
HJCUHX	No blood indicated.	No semen identified.	Saliva indicated.

TABLE 2 - Item 1
Four swabs spotted with 100µl of Saliva.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
HJCUHX	No blood indicated.	No semen identified.	Saliva indicated.
J8DD62	Negative	Negative	Positive
JVLQY7	Negative	Negative	Positive
K28BTQ	Negative	Negative	Positive
K3WX4T	Blood not detected	Semen not detected	Amylase, indicative of saliva
KEA3BP	Negative	Negative	Positive
KETDZ6	Human blood indicated	Negative	Possible presence of saliva
KGHKK2	Negative for blood	Negative for semen	Pos. Human Salivary Amylase
KJN692	Negative	Negative	Positive
L9N8AL	Negative	Negative	Positive
LQLPYG	Negative	Negative	Positive
LRV4H3	Positive	Negative	Possible saliva
LUKB3Z	Negative	Negative	Positive
M3B3FN	Positive for higher primate or ferret blood.	Negative for semen	Positive for α-amylase
M48Q7U	No blood indicated.	No semen identified.	Saliva indicated.
MGAGMT	Negative	Negative	Positive
N4Y7VW		Negative	Positive
NALRPG	Negative	Negative	Indicative
PF8YRL	Blood not detected	Semen not detected	Saliva possibly detected
R7EA6H	Negative	Negative	Positive

TABLE 2 - Item 1
Four swabs spotted with 100µl of Saliva.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
R7YFFN	No blood indicated.	No semen identified.	Saliva indicated
RJZCCV	Human blood indicated	Negative	Possible presence indicated
T9XLJE	Negative	Negative	
TGBTYH	Negative	Negative	Positive
TGTPMB	Negative	Negative	Positive
ULKDJ9	Negative	Negative	Positive
VK47JR	No blood indicated	No semen identified.	Saliva indicated
WAZNX9	No blood present	No semen present	Saliva present
WLAQ8R		Negative	Positive
WQ4H49	No blood was indicated.	No seminal fluid was indicated.	
WQNYJB	Blood not detected	Semen not detected	Amylase, indicative of saliva
XT8BX7	Neg (=) No indication for blood	Neg (=)	Inconclusive
XXLHBR	Human blood indicated	Negative	Negative
Y7WPEE	Blood not detected	Semen not detected	Amylase, indicative of saliva
ZG8GHE	Blood not detected	Semen not detected	Amylase, indicative of saliva
ZHYHGB	Negative	Negative	Positive
ZVKJNN	Positive	Negative	*Positive - see comment page 9[Table 3]
ZW33KZ	Negative	Negative	Positive

TABLE 2 - Item 2
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
2NMLP6	Positive	Negative	Negative
2XTJD9	Positive	Negative	
4F8HWC	Human blood detected**	Semen not detected	Saliva not detected
6DCQBC	Blood indicated.	No semen identified.	No saliva indicated.
6RBUDH	Positive for human blood	Negative	Negative
7A8K9C	Blood indicated.	No semen identified.	No saliva indicated.
7N9F6H	+ positive	- negative	possible
7T4D8A	Human Blood Detected	Semen not detected	Amylase - Various Bodily Fluid
7Y6RN4	Positive	Negative	Negative
8D22XB	Blood indicated	No semen identified	No saliva indicated
8WXUW9	positive. apparent blood present	Negative	Negative
927828	blood indicated	No semen identified	No saliva indicated
AWM2QV	Positive	Negative	Negative
BLJ79Q	Positive	Negative	Negative
BUXV7Y	Positive	Negative	Negative
C74XV9	Positive	Negative	Negative
CT4Y9C	Blood indicated.	No semen identified.	No saliva indicated.
CUWH34	Human blood detected	Semen not detected	Amylase - various body fluids
CWLNMZ	Positive	Negative	Positive
CY79E6	Blood indicated.	No semen identified.	No saliva indicated.

TABLE 2 - Item 2
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
D8Y9K2	Human blood detected	Semen not detected	Amylase not detected
DBGZH6	Blood indicated.	No semen identified.	No saliva indicated.
DBYGLQ	Positive	Negative	Negative
EBXDL8	Positive for a presumptive test for blood	Negative	Negative
ELA9PQ	Human blood indicated	no semen identified	
EXNXE3	Blood indicated	No semen identified	No saliva indicated
F92HC4	Blood present	No semen detected	No saliva detected but unable to see reaction due to supernatant colour, therefore inconclusive.
FB83Y4	Positive	Negative	Negative
FFJD6X	Human blood detected	Semen not detected	Amylase, not indic. of saliva
FGBE6U	Positive	Negative	Negative
FPTPHV	Positive	Negative	Negative
GE84NK	Positive (see comments section)	Negative	Negative
GUWFC9	Positive	Negative	Possible
GZ76ZN	Positive: indicative for human blood	Negative	Negative
H4ZKQM	human blood was indicated	no semen identified	
H7PRBJ	Positive	Negative	Negative
HAPXUW	Human blood detected	Semen not detected	Amylase not detected
HJCUHX	Blood indicated.	No semen identified.	No saliva indicated.
J8DD62	Positive	Negative	Negative

TABLE 2 - Item 2
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
J8DD62	Positive	Negative	Negative
JVLQY7	Positive	Negative	Negative
K28BTQ	Positive	Negative	Negative
K3WX4T	Human blood detected	Semen not detected	Amylase not detected
KEA3BP	Positive	Negative	Negative
KETDZ6	Human blood indicated	Negative	Possible presence of saliva
KGHKK2	Positive for blood	Negative for semen	Negative for Human Salivary Amylase
KJN692	Presumptive Positive	Negative	Negative
L9N8AL	Positive	Negative	Negative
LQLPYG	Positive	Negative	Negative
LRV4H3	Positive	Negative	Negative
LUKB3Z	Positive	Negative	Negative
M3B3FN	Positive for higher primate or ferret blood	Negative for semen	Negative for a-amylase
M48Q7U	Blood indicated.	No semen identified.	No saliva indicated.
MGAGMT	Human blood positive	Negative	Negative
N4Y7VW	Positive	Negative	Negative
NALRPG	Positive	Negative	Negative
PF8YRL	Blood possibly detected	Semen not detected	Saliva not detected
R7EA6H	Positive	Negative	Negative
R7YFFN	Blood indicated	No semen identified	No saliva indicated

TABLE 2 - Item 2
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
RJZCCV	Human blood indicated	Negative	Possible presence indicated
T9XLJE	Positive	Negative	
TGBTYH	Positive	Negative	Negative
TGTPMB	Positive	Negative	Negative
ULKDJ9	Positive	Negative	Negative
VK47JR	Blood indicated	No semen identified	No saliva indicated.
WAZNX9	Human blood present	No semen present	No saliva present
WLAQ8R	Positive	Negative	Negative
WQ4H49	Human blood was indicated.	No semen was identified	
WQNYJB	Human blood detected	Semen not detected	Amylase not detected
XT8BX7	Pos (+) indication for blood	Neg (=)	Neg (=)
XXLHBR	Human blood indicated	Negative	Possible saliva
Y7WPEE	Human blood detected	Semen not detected	Amylase, not indic. of saliva
ZG8GHE	Human blood detected	Semen not detected	Amylase not detected
ZHYHGB	Positive	Negative	Negative
ZVKJNN	Positive	Negative	Negative
ZW33KZ	Positive	Negative	Negative

TABLE 2 - Item 3
Four swabs spotted with 100µl of Semen.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
2NMLP6	Negative	Positive	Negative
2XTJD9	Negative	Positive	
4F8HWC	Blood not detected	Human semen detected	Saliva not detected
6DCQBC	No blood indicated.	Semen identified.	No saliva indicated.
6RBUDH	Negative	Positive	Negative
7A8K9C	No blood indicated.	Semen identified.	No saliva indicated.
7N9F6H		+ positive	
7T4D8A		Semen was detected	Amylase not detected
7Y6RN4	Negative	Positive	Negative
8D22XB	No blood indicated	Semen identified	No saliva indicated
8WXUW9	Negative	Positive: semen present	Inconclusive
927828	No blood indicated	Semen identified	No saliva indicated
AWM2QV	Negative	Positive	Negative
BLJ79Q	Negative	Positive	Negative
BUXV7Y	Negative	Positive	Negative
C74XV9		Positive	
CT4Y9C	No blood indicated.	Semen identified.	No saliva indicated.
CUWH34		Semen detected	Amylase not detected
CWLNMZ	Negative	Positive	Positive
CY79E6	No blood indicated.	Semen identified.	No saliva indicated.
D8Y9K2	Blood not detected	Semen detected	Amylase, not indic. of saliva

TABLE 2 - Item 3
Four swabs spotted with 100µl of Semen.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
D8Y9K2	Blood not detected	Semen detected	Amylase, not indic. of saliva
DBGZH6	No blood indicated.	Semen identified.	No saliva indicated.
DBYGLQ		Positive	
EBXDL8		Positive for semen	Negative
ELA9PQ	no blood indicated	Semen identified	
EXNXE3	No blood indicated	Semen identified	No saliva indicated
F92HC4	No blood detected	Semen present	No saliva detected
FB83Y4	Negative	Positive	Negative
FFJD6X	Blood not detected	Semen detected	Amylase, not indic. of saliva
FGBE6U	Negative	Positive	Negative
FPTPHV	Negative	Positive	Negative
GE84NK	Negative	Positive (see comments section)	Negative
GUWFC9		Positive	
GZ76ZN	Negative	Positive	Negative
H4ZKQM	no blood was indicated	Semen was identified	
H7PRBJ	Negative	Positive	Positive
HAPXUW	Blood not detected	Semen detected	Amylase not detected
HJCUHX	No blood indicated.	Semen identified.	No saliva indicated.
J8DD62		Positive	
JVLQY7		Positive	

TABLE 2 - Item 3
Four swabs spotted with 100µl of Semen.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
K28BTQ	Negative	Positive	Negative
K3WX4T	Blood not detected	Semen detected	Amylase, not indic. of saliva
KEA3BP	Negative	Positive	Negative
KETDZ6		Positive for Sperm	
KGHKK2	Negative for blood	Positive for semen	Negative for Human Salivary Amylase
KJN692	Negative	Positive	Negative
L9N8AL	Negative	Positive	Negative
LQLPYG	Negative	Positive	Negative
LRV4H3		Positive	
LUKB3Z		Positive	
M3B3FN	Presumptive Negative for blood	Positive for Semen	Negative for a-amylase
M48Q7U	No blood indicated.	Semen identified.	No saliva indicated.
MGAGMT	Human blood positive	Positive	Negative
N4Y7VW		Positive	Negative
NALRPG	Negative	Positive	Negative
PF8YRL	Blood not detected	Semen detected	Saliva not detected
R7EA6H	Negative	Positive	Negative
R7YFFN	No blood indicated	Semen identified	No saliva indicated
RJZCCV		Sperm +, possible seminal fluid	
T9XLJE		Positive	

TABLE 2 - Item 3
Four swabs spotted with 100µl of Semen.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
T9XLJE		Positive	
TGBTYH	Negative	Positive	Positive
TGTPMB		Positive	
ULKDJ9	Negative	Positive	Negative
VK47JR	No blood indicated	Semen identified.	No saliva indicated
WAZNX9	No blood present	Semen present	No saliva present
WLAQ8R		Positive	Negative
WQ4H49	No blood was indicated	Semen was identified.	
WQNYJB	Not detected	Semen detected	Amylase not detected
XT8BX7		Pos (+)	
XXLHBR		Positive	
Y7WPEE	Blood not detected	Semen detected	Amylase, not indic. of saliva
ZG8GHE	Blood not detected	Semen detected	Amylase not detected
ZHYHGB	Negative	Positive	Negative
ZVKJNN		*Positive - see comment page 9[Table 3]	
ZW33KZ	Negative	Positive	Positive

TABLE 2 - Item 4
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
2NMLP6	Positive	Negative	Negative
2XTJD9	Positive	Negative	
4F8HWC	Human blood detected**	Semen not detected	Saliva not detected
6DCQBC	Blood indicated.	No semen identified.	No saliva indicated.
6RBUDH	Positive for human blood	Negative	Negative
7A8K9C	Blood indicated.	No semen identified.	No saliva indicated.
7N9F6H	+ positive		
7T4D8A	Human Blood Detected		
7Y6RN4	Positive	Negative	Negative
8D22XB	Blood indicated	No semen identified	No saliva indicated
8WXUW9	positive. apparent blood presen[sic]	Negative	Negative
927828	Blood indicated	No semen identified	No saliva indicated
AWM2QV	Positive	Negative	Negative
BLJ79Q	Positive	Negative	Negative
BUXV7Y	Positive	Negative	Negative
C74XV9	Positive	Negative	
CT4Y9C	Blood indicated.	No semen identified.	No saliva indicated.
CUWH34	Human blood detected		
CWLNMZ	Positive	Negative	Positive
CY79E6	Blood indicated.	No semen identified.	No saliva indicated.

TABLE 2 - Item 4
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
D8Y9K2	Human blood detected	Semen not detected	Amylase not detected
DBGZH6	Blood indicated.	No semen identified.	No saliva indicated.
DBYGLQ	Positive		
EBXDL8	Positive for a presumptive test for blood	Negative	Negative
ELA9PQ	Human blood indicated	no semen identified	
EXNXE3	Blood indicated	No semen identified	No saliva indicated
F92HC4	Blood found	No semen detected	Inconclusive - unable to see reaction due to colour of supernate.
FB83Y4	Positive	Negative	Negative
FFJD6X	Human blood detected	Semen not detected	Amylase, not indic. of saliva
FGBE6U	Positive	Negative	Negative
FPTPHV	Positive	Negative	Negative
GE84NK	Positive	Negative	Negative
GUWFC9	Positive	Positive seminal fluid	
GZ76ZN	Positive	Negative	Negative
H4ZKQM	human blood was indicated	no semen was identified	
H7PRBJ	Positive	Negative	Negative
HAPXUW	Human blood detected	Semen not detected	Amylase not detected
HJCUHX	Blood indicated.	No semen identified.	No saliva indicated.
J8DD62	Positive	Negative	
JVLQY7	Positive	Negative	

TABLE 2 - Item 4
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
JVLQY7	Positive	Negative	
K28BTQ	Positive	Negative	Negative
K3WX4T	Human blood detected	Semen not detected	Amylase not detected
KEA3BP	Positive	Negative	Negative
KETDZ6	Human blood indicated	Seminal fluid indicated	
KGHKK2	Positive for blood	Negative for semen	Negative for Human Salivary Amylase
KJN692	Presumptive Positive	Negative	Negative
L9N8AL	Positive	Negative	Negative
LQLPYG	Positive	Negative	Negative
LRV4H3	Positive	Negative	
LUKB3Z	Positive	Negative	
M3B3FN	Positive for higher primate or ferret blood	Negative for semen	Negative for α-amylase
M48Q7U	Blood indicated.	No semen identified.	No saliva indicated.
MGAGMT	Human blood positive	Negative	Negative
N4Y7VW	Positive	Negative	Negative
NALRPG	Positive	Negative	Negative
PF8YRL	Blood possibly detected	Semen not detected	Saliva not detected
R7EA6H	Positive	Negative	Negative
R7YFFN	Blood indicated	No semen identified	No saliva indicated
RJZCCV	Human blood indicated	Negative	

TABLE 2 - Item 4
Four swabs spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
T9XLJE	Positive	Indicated only	
TGBTYH	Positive	Negative	Negative
TGTPMB	Positive		
ULKDJ9	Positive	Negative	Negative
VK47JR	Blood indicated	No semen identified.	No saliva indicated
WAZNX9	Human blood present	No semen present	No saliva present
WLAQ8R	Positive	Negative	Negative
WQ4H49	Human blood was indicated.	No semen was identified.	
WQNYJB	Human blood detected	Semen not detected	Amylase not detected
XT8BX7	Pos (+) indications for blood		
XXLHBR	Human blood indicated	Negative	
Y7WPEE	Human blood detected	Semen not detected	Amylase, not indic. of saliva
ZG8GHE	Human blood detected	Semen not detected	Amylase not detected
ZHYHGB	Positive	Negative	Negative
ZVKJNN	Positive		
ZW33KZ	Positive	Negative	Negative

TABLE 2 - Item 5
One glove spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
2NMLP6	Positive	Negative	Negative
2XTJD9	Positive	Negative	
4F8HWC	Human blood detected**	Semen not detected	Saliva not detected
6DCQBC	Blood indicated.	No semen identified.	No saliva indicated.
6RBUDH	Positive for human blood	Negative	
7A8K9C	Blood indicated.	No semen identified.	No saliva indicated.
7N9F6H	+ positive	- negative	
7T4D8A	Human blood detected		Amylase - various body fluids
7Y6RN4	Positive	Negative	Negative
8D22XB	Blood indicated	No semen identified	No saliva indicated
8WXUW9	positive. apparent blood present		N/A no conclusions with ALS
927828	Blood indicated	No semen identified	No saliva indicated
AWM2QV	Positive	Negative	Negative
BLJ79Q	Positive	Negative	
BUXV7Y	Positive	Negative	Negative
C74XV9	Positive		
CT4Y9C	Blood indicated.	No semen identified.	No saliva indicated.
CUWH34	Human blood detected		Amylase - various body fluids
CWLNMZ	Positive	Negative	Positive
CY79E6	Blood indicated.	No semen identified.	No saliva indicated.

TABLE 2 - Item 5
One glove spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
D8Y9K2	Human blood detected		Amylase not detected
DBGZH6	Blood indicated.	No semen identified.	No saliva indicated.
DBYGLQ	Positive	Negative	
EBXDL8	Positive for a presumptive test for blood	Negative	Negative
ELA9PQ	Human blood indicated	no semen identified	
EXNXE3	Blood indicated	No semen identified	No saliva indicated
F92HC4	Blood present on inside & outside		
FB83Y4	Positive	Negative	Negative
FFJD6X	Human blood detected		Amylase not detected
FGBE6U	Positive	Negative	Negative
FPTPHV	Positive	Negative	Negative
GE84NK	Positive	Negative	Negative
GUWFC9	Positive		
GZ76ZN	Positive	Negative	Negative
H4ZKQM	human blood was indicated	no semen was identified	
H7PRBJ	Positive	Negative	Negative
HAPXUW	Human blood detected		Amylase not detected
HJCUHX	Blood indicated.	No semen identified.	No saliva indicated.
J8DD62	Positive		
JVLQY7	Positive		

TABLE 2 - Item 5
One glove spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
K28BTQ	Positive	Negative	Negative
K3WX4T	Human blood detected		Amylase not detected
KEA3BP	Positive	Negative	Negative
KETDZ6	Human blood indicated		
KGHKK2	Positive for blood	Negative for semen	Neg. for Human Salivary Amylase
KJN692	Presumptive Positive	Negative	Negative
L9N8AL	Positive	Negative	Negative
LQLPYG	Positive	Negative	Negative
LRV4H3	Positive		
LUKB3Z	Positive		
M3B3FN	Positive for higher primate or ferret blood	Presumptive negative for semen	
M48Q7U	Blood indicated.	No semen identified.	No saliva indicated.
MGAGMT	Human blood positive	Negative	Negative
N4Y7VW	Positive	Negative	Negative
NALRPG	Positive		
PF8YRL	Blood possibly detected	Semen not detected	Saliva not detected
R7EA6H	Positive	Negative	Negative
R7YFFN	Blood indicated	No semen identified	No saliva indicated
RJZCCV	Human blood indicated	No stains visualized using ALS	
T9XLJE	Positive		

TABLE 2 - Item 5
One glove spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
T9XLJE	Positive		
TGBTYH	Positive	Negative	
TGTPMB	Positive		
ULKDJ9	Positive	Negative	Negative
VK47JR	Blood indicated.	No semen identified.	No saliva indicated.
WAZNX9	Human blood present	No semen present	No saliva present
WLAQ8R	Positive	Negative	Negative
WQ4H49	Human blood was indicated.	No semen was identified.	
WQNYJB	Human blood detected		Amylase not detected
XT8BX7	Pos (+) indication for blood		
XXLHBR			
Y7WPEE	Human blood detected		Amylase not detected
ZG8GHE	Human blood detected		Amylase not detected
ZHYHGB	Positive	Negative	Negative
ZVKJNN	Positive	*Negative - see comment page 9[Table 3]	*Negative - see comment page 9[Table 3]
ZW33KZ	Positive	Negative	Negative

TABLE 2 - Item 6
One white tank top spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
2NMLP6	Positive	Negative	Negative
2XTJD9	Positive	Positive	
4F8HWC	Human blood detected**	Semen not detected	Saliva not detected
6DCQBC	Blood indicated.	No semen identified.	No saliva indicated.
6RBUDH	Positive for human blood	Negative	
7A8K9C	Blood indicated.	No semen identified.	No saliva indicated.
7N9F6H	+ positive	- negative	
7T4D8A	Human blood detected		Amylase - various body fluids
7Y6RN4	Positive	Negative	Negative
8D22XB	Blood indicated	No semen identified	No saliva indicated
8WXUW9	Positive	Negative	N/A. entire item not tested.
927828	Blood indicated	No semen identified	No saliva indicated
AWM2QV	Positive	Negative	Negative
BLJ79Q	Positive	Negative	
BUXV7Y	Positive	Negative	Negative
C74XV9	Positive	Negative	
CT4Y9C	Blood indicated.	No semen identified.	No saliva indicated.
CUWH34	Human blood detected	Semen not detected	Amylase - various body fluids
CWLNMZ	Positive	Negative	Positive
CY79E6	Blood indicated.	No semen identified.	No saliva indicated.
D8Y9K2	Human blood detected		Amylase not detected

TABLE 2 - Item 6
One white tank top spotted with 100 μ l of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
D8Y9K2	Human blood detected		Amylase not detected
DBGZH6	Blood indicated.	No semen identified.	No saliva indicated.
DBYGLQ	Positive	Negative	
EBXDL8	Positive for a presumptive test for blood	Negative	Negative
ELA9PQ	Human blood indicated	no semen identified	
EXNXE3	Blood indicated	No semen identified	No saliva indicated
F92HC4	Blood present	No semen detected	
FB83Y4	Positive	Negative	Negative
FFJD6X	Human blood detected		Amylase not detected
FGBE6U	Positive	Negative	Negative
FPTPHV	Positive	Negative	Negative
GE84NK	Positive	Negative	Negative
GUWFC9	Positive		
GZ76ZN	Positive	Negative	Negative
H4ZKQM	human blood was indicated	no semen was identified	
H7PRBJ	Positive	Negative	Negative
HAPXUW	Human blood detected		Amylase not detected
HJCUHX	Blood indicated.	No semen identified.	No saliva indicated.
J8DD62	Positive	Negative	
JVLQY7	Positive		

TABLE 2 - Item 6
One white tank top spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
K28BTQ	Positive	Negative	Negative
K3WX4T	Human blood detected		Amylase not detected
KEA3BP	Positive	Negative	Negative
KETDZ6			
KGHKK2	Positive for blood	Negative for semen	Neg. for Human Salivary Amylase
KJN692	Presumptive Positive	Negative	Negative
L9N8AL	Positive	Negative	Negative
LQLPYG	Positive	Negative	Negative
LRV4H3	Positive		
LUKB3Z	Positive	Negative	
M3B3FN	Positive for higher primate or ferret blood	Presumptive negative for semen	
M48Q7U	Blood indicated.	No semen identified.	No saliva indicated.
MGAGMT	Human blood positive	Negative	Negative
N4Y7VW	Positive	Negative	Negative
NALRPG	Positive		
PF8YRL	Blood possibly detected	Semen not detected	Saliva not detected
R7EA6H	Positive	Negative	Negative
R7YFFN	Blood indicated	No semen identified	No saliva indicated
RJZCCV	Human blood indicated	No stains visualized using ALS	
T9XLJE	Positive		

TABLE 2 - Item 6
One white tank top spotted with 100µl of Blood.

Response Summary			
Webcode	Blood Conclusion	Semen Conclusion	Saliva Conclusion
T9XLJE	Positive		
TGBTYH	Positive	Negative	
TGTPMB	Positive		
ULKDJ9	Positive	Negative	Negative
VK47JR	Blood indicated.	No semen identified.	No saliva indicated
WAZNX9	Human blood present	No semen present	No saliva present
WLAQ8R	Positive	Negative	Negative
WQ4H49	Human blood was indicated.	No semen was identified.	
WQNYJB	Human blood detected		Amylase not detected
XT8BX7	Pos (+) indication for Blood		
XXLHBR	Human Blood Indicated		
Y7WPEE	Human blood detected		Amylase not detected
ZG8GHE	Human blood detected		Amylase not detected
ZHYHGB	Positive	Negative	Negative
ZVKJNN	Positive	*Negative - see comment page 9[Table 3]	*Negative - see comment page 9[Table 3]
ZW33KZ	Positive	Negative	Negative

Additional Comments

TABLE 3

WebCode	Additional Comments
4F8HWC	** Reports state 'human blood detected' if Haematrace[sic] positive, however, a glossary is attached to the report explaining the possible cross reactivity with the blood of apes and domestic ferrets. [No glossary was provided by participant.] *Positive results with the Phadebas amylase test are reported as the 'possible presence of saliva'. If nucleated squamous epithelial cells are found on a slide prepared from the positive area or swab, these two results in conjunction may be reported in a stronger manner, depending on the rapidity and strength of the Phadebas test, number of cells found, case information etc.
7N9F6H	Item #1, swabs with red paint were tested with LMG and gave negative results. In actual case work, analyst would not have progressed to Hematrace. Our report requires a statement to indicate that test results were out of the ordinary: "confirmatory test was positive for indication of human blood; however, preliminary testing was negative."
8WXUW9	All scenario details (which includes the collection location of swabs and items) should be written on the first page of the document, under the heading 'Scenario'. This is especially important when a comment such as 'taken from the locations described in the above scenario' is written next to each of the samples in the 'Items submitted' list. The additional information (swab of the victims hands, swab from the inner thigh, etc.) which was scattered on each of the following pages of the 10 page document is easily missed. There really isn't enough pertinent details included in the scenario. In the real world, if there is missing information, we can contact the submitting agency and request this information as it is pivotal to how items are examined. For example, is the victim bleeding? and how was the victim murdered? strangled/suffocated only, blunt/sharp force trauma, shooting? This affects what body fluids are looked for, on what items, and to what degree. With shooting and blunt force trauma, you may only find a small <1mm stain on a garment and thus these cases may require a much more in-depth analysis of items. Under regular case circumstances, I wouldn't screen a suspects tank top or gloves for saliva or semen, unless there are some specific circumstances which warrants it (e.g. victim was suffocated), but given this is a CTS and a body fluid identification test, some of these were conducted anyway. Trace DNA would/may be collected from some of these items (e.g. glove) instead of wasting potential sample during testing.
CWLNMZ	RSID - (Rapid Stain Identification Kit) manufactured by Independent forensics DNA Testing and Technologies, is the kit used in this laboratory for confirmatory testing of blood, saliva, and semen.
EBXDL8	The six items were examined for the presence of blood, semen and saliva. Only one body fluid was identified on each of the six items examined. At [Laboratory] there is no confirmatory test available for the presence of blood or saliva. As a result our statements state: For saliva: the item was examined for the presence of amylase, a constituent of saliva. For blood: item was examined for blood and there were areas present on the item which gave a positive reaction to a presumptive test for blood. For item 3 the presence of semen was confirmed microscopically. Further work: The samples could be submitted for DNA analysis along with reference[sic] samples from the victim and the suspect to attribute the body fluids to a person.
EXNXE3	Please note that all tests were performed on all items (regardless of scenario) per internal agency instructions. Semen reporting is based on all tests (i.e., item #1 had +1 AP result but "no semen identified" report wording is based on AP [+1], PSA [neg], and Micro [neg]). ALS was used to screen Items #5 and #6 for fluorescing stains. Any fluorescing stains were then tested for semen/saliva using the tests as described on the answer sheet. Please reference notes for further information if needed. [Reference notes were not provided by this participant.]

TABLE 3

WebCode	Additional Comments
GE84NK	Presumptive chemical tests for the presence of blood is done with ortho-tolidine. A positive presumptive chemical test result does not confirm the presence of blood. Testing for the presence of human blood is done with ABACard HeamTrace[sic] test kit. This kit is known to cross-react with some species of the Mustelidae family and with higher primates. Presumptive chemical testing for the presence of seminal fluid/semen is done by the Modified Walker Acid Phosphatase test. A positive presumptive chemical test does not confirm the presence of seminal fluid/semen. Testing for the presence of PSA is done with the ABACard P30 test kit. This kit may yield positive results with male urine. Sperm cell identification is done by microscopic examination of stained (nuclear fast red/picroindigocarmine) slides. Testing for the presence of human salivary alpha-amylase (an enzyme found in human saliva) is done with the Rapid Stain Identification-Saliva test kit (RSID-Saliva). This test is known to yield weakly positive results with breast milk and fecal samples.
GUWFC9	Item 3: Blood test not performed because no red brown staining visualized. Saliva test not performed to conserve sample. Item 4: Presumptive test for seminal fluid not performed due to color of stain. Saliva test not performed to conserve sample. Item 5: Interior and exterior of glove swabbed for the possible presence of epithelial cells. Item 6: Interior collar swabbed for the possible presence of epithelial cells.
GZ76ZN	Re Clothing found in the suspect's: Blood was found on the glove and on the top. During the course of normal casework this staining would have been submitted for DNA profiling. The Haematracesic test would only be carried out depending on case circumstances and the DNA result. The Haematracesic test was not carried out on (item 4) as the swabs were taken from the victims hand. If the items of clothing belonged to the suspect these items would not have been screened for the presence of semen as the semen source may be the suspect's own. Even though the items of clothing were recovered from the suspect's home, I am not sure who owns them. To establish the wearer, the inside of both items were minitaped for the presence of epithelial cells. Minitaping was carried out on the outside of the glove in the event that if epithelial cells were transferred from someone else/touch DNA (e.g. victim) other than the wearer. Fibre lifting would be considered depending on the DNA result. Fibre Lifting was carried out on the top. Def: Minitape lifting is a process whereby adhesive strips remove epithelial (skin cells) that may be present on surfaces of items.
H7PRBJ	A positive presumptive chemical test result does not confirm the presence of blood. Presumptive chemical testing for the presence of blood is done with ortho-Tolidine. Testing for the presence of human blood is done with ABACard HemaTrace test kit. This kit is known to cross react with some species of the Mustelidae family and with higher primates. Presumptive chemical testing for the presence of seminal fluid/semen is done by the Modified Walker Acid Phosphatase test. A positive presumptive chemical test does not confirm the presence of seminal fluid/semen. Testing for the presence of PSA is done with the ABACard P30 test kit. This kit may yield positive results with male urine. Sperm cell identification is done by microscopic examination of stained (nuclear fast red/picroindigocarmine) slides. Testing for the presence of human salivary alpha-amylase (an enzyme found in human saliva) is done with the Rapid Stain Identification-Saliva test kit (RSID-Saliva). This kit is known to yield weakly positive results with breastmilk and fecal samples.
KETDZ6	Item #1: Epithelial cells were microscopically visualized on a slide made from the sample. Item #2: The possible presence of epithelial cells was noted. Item #3: Per laboratory protocol sample was reported to be positive for sperm. Presence of seminal fluid not confirmed due to microscopic results. Item #5: Swabs were taken from the interior cuff of the glove for the possible presence of epithelial cells. A suspected pink fiber was also recovered from the interior "ring finger" area of the glove and stored for future testing if necessary. Item #6: Swabs were taken from the interior collar and waistband area of the tank top for the possible presence of epithelial cells.
KGHKK2	Note: Actual conclusion for blood is: "Presumptive phenolphthalein testing was positive for the presence of blood".

TABLE 3

WebCode	Additional Comments
KJN692	The laboratory does not conduct any tests for the confirmation of blood. Report conclusions on a "dried red stain" are written as follows: Presumptive phenolphthalein testing for the presence of blood was positive (or negative) on the... Item #6 also contained two off white / yellow stains each approx 1cm in diameter on the upper front of the shirt. The stains were visible under normal light. The stains were negative ALS, KM and AP. It is assumed that these stains were not part of the actual test and inadvertently ended up on the item. The characteristic appearance of blood is not an actual test, it [sic] only based on the initial visual examination of an item.
LQLPYG	Presumptive chemical testing for the presence of blood is done with ortho-Tolidine. A positive presumptive chemical test result does not confirm the presence of blood. Testing for the presence of human blood is done with the ABACard HemaTrace test kit. This kit is known to cross react with some species of the Mustelidae family and with higher primates. Presumptive chemical testing for the presence of seminal fluid / semen is done by the Modified Walker Acid Phosphatase test. A positive presumptive chemical test result does not confirm the presence of seminal fluid / semen. Testing for the presence of PSA is done with the ABACard p30 test kit. This kit may yield positive results with male urine. Sperm cell identification is done by microscopic examination of stained (nuclear fast red / picroindigocarmin) slides. Testing for the presence of human salivary alpha-amylase (an enzyme found in human saliva) is done with the Rapid Stain Identification-Saliva test kit (RSID-Saliva). This kit is known to yield weakly positive results with breastmilk and fecal samples.
M48Q7U	The control swabs were also tested with the following tests with the following results: KM - Neg; AP - Neg; PSA - Neg; Micro - KPIC - Neg; Phadebas - Neg. Overall conclusions: Blood - No blood indicated. Semen - No semen identified. Saliva - No saliva indicated. Please note that all stains from Items 1-6 were tested for each body fluid per direction of the evidence control center manager, regardless of the scenario given.
NALRPG	ALS was used as a screening technique. As such, it is neither a presumptive nor confirmatory test. Item 2 and control swabs were tested for the presence of the following animal proteins: bovine, cat, chicken, deer, dog, swine, horse, sheep, and rabbit. All results were negative.
RJZCCV	Item 4: Presumptive screening (AP) for seminal fluid was not performed due to the staining on the swabs. Item 5: The item was swabbed for the possible presence of epithelial cells. Item 6: The item was swabbed for the possible presence of epithelial cells.
TGTPMB	ALS used to view Item 5
ULKDJ9	Presumptive chemical testing for the presence of blood is done with ortho-Tolidine. A positive presumptive chemical test result does not confirm the presence of blood. Testing for the presence of human blood is done with the ABACard HemaTrace test kit. This kit is known to cross react with some species of the Mustelidae family and with higher primates. Presumptive chemical testing for the presence of seminal fluid/semen is done by the Modified Walker Acid Phosphatase test. A positive presumptive chemical test does not confirm the presence of seminal fluid/semen. Testing for the presence of PSA is done with the ABACard p30 test kit. This kit may yield positive results with male urine. Testing for the presence of human salivary alpha-amylase (an enzyme found in human saliva) is done with Rapid Stain Identification-Saliva test kit (RSID-Saliva). This test is known to yield weakly positive results with breast milk and fecal samples.
WLAQ8R	Only two control swaps[sic] for 6 items are not enough! The test is too easy. There should be more than one body fluid on a swab for example.
WQ4H49	No hairs were observed on the submitted items.
XXLHBR	Items 5 and 6: A collection for epithelial cells was completed. No testing was performed to verify the presence of cells, the sample would be sent to DNA for further analysis. Item 5: The glove was screened for the possible presence of saliva and seminal fluid using an ALS.

TABLE 3

WebCode	Additional Comments
ZVKJNN	<p>Item 1: (Saliva Tests) The conclusion according to [Laboratory] reporting results for positive starch iodine test for saliva would be possible saliva since we are testing for alpha amylase. However for the Summary Report as requested by CTS the result for saliva was recorded as positive. Item 3: (Semen Test) The conclusion according to [Laboratory] reporting results for positive AP test without p30 confirmation is the possible presence of seminal fluid. The identification of sperm is reported as positive for sperm. However for the Summary Report as requested by CTS the result for semen was recorded as positive. Items 5 & 6: Per procedure at [Laboratory] these items were swabbed on the interior for the possible presence of epithelial cells. No microslide was made, items are forwarded to DNA for further testing. Although other body fluids can fluoresce using an ALS, semen and saliva are the most common body fluids that would fluoresce that [Laboratory] would test for.</p>
ZW33KZ	<p>Semen was identified on Item 3. Tests for salivary alpha-amylase, a component of saliva, were weakly positive on this item. The results of these tests indicate a mixture of semen and saliva on Item 3. Testing for the presence of human salivary alpha-amylase, an enzyme found in human saliva, is done with Rapid Stain Identification-Saliva test kit (RSID-Saliva). This test is known to yield weakly positive results with breast milk and fecal samples. Testing for the presence of PSA is done with ABACard p30. This kit may yield positive results with male urine. Sperm cell identification is done by microscopic examination of stained (nuclear fast red/picroindigocarmine) slides. Testing for the presence of human blood is done with the ABACard HemaTrace test kit. This kit is known to cross react with some species of the Mustelidae family and with higher primates.</p>

Appendix: Data Sheet

Collaborative Testing Services ~ Forensic Testing Program

Test No. 11-579: Body Fluid Identification Test

DATA MUST BE RECEIVED BY November 7, 2011 TO BE INCLUDED IN THE REPORT

Participant Code:

WebCode:

Accreditation Release Statement

CTS submits external proficiency test data directly to ASCLD/LAB and FQS-International. Please select one of the following statements to ensure your data is handled appropriately.

- This participant's data is intended for submission to ASCLD/LAB and/or FQS-International. (Accreditation Release section on the last page must be completed and submitted.)
- This participant's data is **NOT** intended for submission to ASCLD/LAB or FQS-International.

Scenario (this scenario applies only to Items 3 - 6):

Police are investigating the rape and murder of a woman found dead in her apartment. The landlord stated that there was a man loitering outside her apartment earlier that morning. A description of the man led investigators to the suspect. Investigators have submitted swabs gathered from the body of the victim and clothing items recovered from the suspect's home.

Items Submitted (Sample Pack B9):

- Items 1-2 Swabs (4 each) - scenario not applicable
- Items 3-4 Swabs (4 each) - taken from locations described in the above scenario
- Items 5-6 Clothing items - taken from locations described in the above scenario
- 1 package of unused swabs

Please return all pages of this data sheet.

Page 1 of 10

Instructions:

Please provide the results of any screening tests for body fluids performed on the items submitted on the appropriate page for each Item. List the tests performed, the results of each test, and indicate whether the test was considered presumptive or confirmatory. Item descriptions (as applicable) are provided at the top of each Item page.

For the purposes of this test, all swabs within an Item are to be treated as duplicate swabs. A package of unused swabs is also provided for control purposes.

Please use the abbreviations listed in this response key to fill in the Screening Test tables. This is not an all inclusive list of tests, and should not be used to determine what tests should be performed. TESTS NOT ON THIS LIST MAY BE USED FOR SCREENING.

Test	Abbreviation	Test	Abbreviation
Acid Phosphatase	AP	Phenolphthalein-Tetramethyl benzidine	PTMB
Kastle Meyer	KM	Prostate Specific Antigen	PSA
Leucomalachite Green	LMG	Alternate Light Source	ALS
Microscopic	Micro	Tetramethyl benzidine	TMB
Ortho-tolidine	O-tol		

Please list the tests performed for each type of body fluid screened for and on the corresponding line, indicate the response with a check in the appropriate box.

CTS recognizes that the "Test Results" wording may differ from your laboratory's terminology, however, to facilitate the Summary Reports, please use these options as best you can.

Test(s) Performed	Test Results			Presumptive / Confirmatory	
	Pos.	Neg.	Inc.		
Hematrace	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
KM	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

****** Indicate Presumptive or Confirmatory for each test.

What is the overall conclusion for blood?
Positive

** The questions in this section are for informational purposes only. These results will not be submitted to Accrediting bodies.

Item 1: Swabs marked with red paint

Scenario: Scenario is not applicable for this Item. Please examine this item for the body fluids that you typically screen for in casework.

<p>Blood Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for blood?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Semen Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for semen?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Saliva Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for saliva?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Other Tests</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:15%;">Body Fluid</th> <th style="width:20%;">Test(s) Performed</th> <th colspan="3" style="text-align: center;">Test Results</th> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;">Pos.</td> <td style="text-align: center;">Neg.</td> <td style="text-align: center;">Inc.</td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Body Fluid	Test(s) Performed	Test Results			_____	_____	Pos.	Neg.	Inc.	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Body Fluid	Test(s) Performed	Test Results																																					
_____	_____	Pos.	Neg.	Inc.																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
<input type="checkbox"/>	<input type="checkbox"/>																																						
<input type="checkbox"/>	<input type="checkbox"/>																																						
<input type="checkbox"/>	<input type="checkbox"/>																																						
<input type="checkbox"/>	<input type="checkbox"/>																																						

** The information reported in this section will not be submitted to Accrediting bodies.

Item 2: Swabs marked with yellow paint

Scenario: Scenario is not applicable for this Item. Please examine this item for the body fluids that you typically screen for in casework.

<p>Blood Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <p>Pos. Neg. Inc.</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p>Presumptive / Confirmatory</p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p>What is the overall conclusion for blood?</p> <p>_____</p>
---	---	--

<p>Semen Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <p>Pos. Neg. Inc.</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p>Presumptive / Confirmatory</p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p>What is the overall conclusion for semen?</p> <p>_____</p>
---	---	--

<p>Saliva Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <p>Pos. Neg. Inc.</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/></p>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p>Presumptive / Confirmatory</p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p>What is the overall conclusion for saliva?</p> <p>_____</p>
--	---	---

<p>Other Tests</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:15%;">Body Fluid</th> <th style="width:35%;">Test(s) Performed</th> <th colspan="3" style="text-align: center;">Test Results</th> </tr> <tr> <td></td> <td></td> <th style="text-align: center;">Pos.</th> <th style="text-align: center;">Neg.</th> <th style="text-align: center;">Inc.</th> </tr> <tr> <td>_____</td> <td>_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>_____</td> <td>_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>_____</td> <td>_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>_____</td> <td>_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>					Body Fluid	Test(s) Performed	Test Results					Pos.	Neg.	Inc.	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p>Presumptive / Confirmatory</p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p> <p><input type="checkbox"/> <input type="checkbox"/></p>	
Body Fluid	Test(s) Performed	Test Results																																		
		Pos.	Neg.	Inc.																																
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																

** The information reported in this section will not be submitted to Accrediting bodies.

Item 3: Swabs marked with green paint

Scenario: The swabs marked with green paint were taken from inner thigh on victim.

<p>Blood Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for blood?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Semen Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for semen?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Saliva Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test Results</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for saliva?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Other Tests</p> <table style="width:100%; border-collapse: collapse;"> <tr> <th style="width:15%;">Body Fluid</th> <th style="width:20%;">Test(s) Performed</th> <th colspan="3" style="text-align: center;">Test Results</th> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;">Pos.</td> <td style="text-align: center;">Neg.</td> <td style="text-align: center;">Inc.</td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;">_____</td> <td style="text-align: center;">_____</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Body Fluid	Test(s) Performed	Test Results			_____	_____	Pos.	Neg.	Inc.	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; border-collapse: collapse;"> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Body Fluid	Test(s) Performed	Test Results																																					
_____	_____	Pos.	Neg.	Inc.																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
_____	_____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																			
<input type="checkbox"/>	<input type="checkbox"/>																																						
<input type="checkbox"/>	<input type="checkbox"/>																																						
<input type="checkbox"/>	<input type="checkbox"/>																																						
<input type="checkbox"/>	<input type="checkbox"/>																																						

** The information reported in this section will not be submitted to Accrediting bodies.

Item 4: Swabs marked with blue paint

Scenario: The swabs marked with blue paint were taken from the victim's hand.

<p>Blood Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for blood?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Semen Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for semen?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Saliva Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for saliva?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Other Tests</p> <p>Body Fluid</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									

** The information reported in this section will not be submitted to Accrediting bodies.

Item 5: Glove

Scenario: The item of evidence is a glove taken from the suspect's home.

<p>Blood Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for blood?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Semen Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for semen?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Saliva Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for saliva?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Other Tests</p> <p>Body Fluid</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									

** The information reported in this section will not be submitted to Accrediting bodies.

Item 6: White tank top

Scenario: The item of evidence is a white tank top taken from the suspect's home.

<p>Blood Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for blood?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Semen Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for semen?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Saliva Tests</p> <p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table> <p>What is the overall conclusion for saliva?</p> <p>_____</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																							
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																								

<p>Other Tests</p> <p>Body Fluid</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p>Test(s) Performed</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>	<p style="text-align: center;">Test Results</p> <table style="width:100%; text-align: center;"> <tr> <th style="width:33%;">Pos.</th> <th style="width:33%;">Neg.</th> <th style="width:33%;">Inc.</th> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	Pos.	Neg.	Inc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<p>** Indicate Presumptive or Confirmatory for each test.</p> <p style="text-align: center;">Presumptive / Confirmatory</p> <table style="width:100%; text-align: center;"> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </table>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pos.	Neg.	Inc.																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																								
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									
<input type="checkbox"/>	<input type="checkbox"/>																									

** The information reported in this section will not be submitted to Accrediting bodies.

RELEASE OF DATA TO ACCREDITATION BODIES

The following Accreditation Releases will apply only to:

Participant Code:

WebCode:

for Test No. **11-579: Body Fluid Identification**

This release page must be completed and received by **November 7, 2011** to have this participant's submitted data included in the reports forwarded to the respective Accreditation Bodies.

ASCLD/LAB RELEASE

If your lab has been accredited by ASCLD/LAB and you are submitting this data as part of their external proficiency test requirements, have the laboratory's designated individual complete the following.

The information below must be completed in its entirety for the results to be submitted to ASCLD/LAB.

ASCLD/LAB Legacy Certificate No. _____ ASCLD/LAB International Certificate No. _____

Signature _____ Date _____

Laboratory Name _____

Location (City/State) _____

FQS-INTERNATIONAL RELEASE

If your laboratory maintains its accreditation through FQS-International, please complete the following form in its entirety to have your results forwarded.

FQS-International Certificate No. _____

Signature and Title: _____ Date _____

Laboratory Name _____

Location (City/State) _____

Return Instructions

Accreditation Release

Please submit the completed Accreditation Release at the same time as your full data sheet. See Data Sheet Return Instructions on the previous page.

*Questions? Contact us 8 am-4:30 pm EST
Telephone: +1-571-434-1925
email: forensics@cts-interlab.com*